

Product Data Sheet



Product Name *3-D Life* SG-PVA

Catalog Number M81-3




Description *3-D Life* SG-PVA is a component of the *3-D Life* Hydrogel system. It contains polyvinyl alcohol functionalized with thiol-reactive groups (SG-PVA). It can be used together with PEG-Link (Cat. No. L50-1 or L50-3) or CD-Link (Cat. No. L60-1 or L60-3) to set up slow gelling, cell-compatible hydrogels.

10 x CB (pH 7.2) buffer allows control of pH and osmotic conditions during gel formation in the presence of cells.

For more information and instructions, please consult the General Protocol GP-2 "Preparation of *3-D Life* Slow Gelling Hydrogels" and the *3-D Life* Hydrogels User Guide on www.cellendes.com.

Quantity Allows formation of up to 6 ml *3-D Life* Hydrogel depending on the stiffness of the gel, if used together with thiol-functionalized polymers (*3-D Life* PEG-Link or CD-Link).

Components

Material	Quantity	Concentration of reactive groups	Storage
 SG-PVA	3x 170 µl	30 mmol/L	Short term (≤2 months): 4°C; Long Term: -80°C, avoid frequent freeze/thaw cycles
 10 x CB (pH 7.2)	600 µl	n.a.	Short term (≤2 months): 4°C. Long term: -20°C to -80°C.
 Water	4x 1500 µl	n.a.	Room temperature or lower

All materials are filter-sterilized.

Intended for research use only. Not for use in human therapeutic or diagnostic applications.

Distributed by:
Ilex Life Sciences LLC
Tel: (828) 531-9949
Email: info@ilexlife.com
Web: <https://ilexlife.com/>

