



NRC International GmbH



O₂ ANALYSER CATALOG



Your trusted source for accurate oxygen analysis - Nrc International



Discover the cutting-edge NRC Oxygen Analyzers, meticulously crafted with German precision and manufactured to the highest industry standards. At NRC International, we strive to enhance your diving journey by providing unmatched accuracy, reliability, and safety. before.

NRC International Team

Contents

Oxy Spy Nitrox Analyser	04
Pick O ₂ Nitrox Analyser	05
Overview 2 NRC Analysers	06
Sale and delivery	07




Important Safety Warning:

Proper training and utilization of the Nitrox Analyzer are crucial for ensuring safe Nitrox diving practices. At NRC International, we emphasize the importance of receiving adequate training and following recommended maintenance protocols. We cannot assume liability for accidents arising from improper use, failure to adhere to maintenance guidelines, or incorrect application. Your safety is our utmost priority.

OXY SPY ANALYSER

Self-servicing Oxygen analyser



1 ONE DECIMAL DISPLAY

■ Specially developed to accommodate the precision needed for Technical Diving. (Imperial or metric)

2 INTEGRATED MOD CALCULATOR

■ Partial pressure (PPO2) adjustable from 1.2 to 1.8 Bar

3 SELF SERVICING O2 CELL AND BATTERY

4 AUTOMATIC CALIBRATION AND POWER OFF

Specifications

WEIGHT	530 g
WIDTH	13 cm
LENGTH	29,5 cm
HEIGHT	11 cm

How to use it



INCLUDED WITH THE OXY SPY ANALYSER



Storage box



1 Oxygen sensor d-05



1 Battery cr2430




1 Servicing tool



SCAN FOR SERVICE VIDEO

PICK O₂ ANALYSER

Single use Oxygen analyser



1 SMALL AND LIGHTWEIGHT

- Perfect for travelling or busy dive shops

2 2 DIGITS DIGITAL DISPLAY

3 SIMPLE CALIBRATION BUTTON

4 INTERCHANGEABLE EXTRACTION PIECES

- See below accessories

Specifications

WEIGHT	80 g
WIDTH	4,8 cm
LENGHT	5,9 cm
HEIGHT	7,4 cm

How to use it



RECYCLING PROGRAM

The Pic O₂ Analyzer is a single-use device, once the cell expires, the analyzer becomes non-functional. However, we're committed to reducing our environmental impact through a recycling program. By returning your used product, you contribute to our planet-saving efforts. In return, you'll receive a discount on your next Pic O₂ Analyzer. Contact us for more information.



BCD ADAPTOR

ACCESSORIES



TANK ADAPTOR

OVERVIEW NRC ANALYSERS



Features	Pick O2 Analyser	Oxy Spy Analyser
Oxygen / Nitrox Analyser	✓	✓
Replaceable O2 Cell	✗	✓
Recycling program	✓	✗
Medical Grade O2 Cell	✓	✓
Accessories	✓	✗
Display Digits	Integer Display No Decimal Display	One Decimal Display

Our commitment during product development was to offer the perfect analyzer for every user. At the core of our analyzers, such as OxySpy or Pic O₂, lie the exceptional oxygen sensors. These sensors are expertly installed and precisely calibrated during manufacturing, a feature shared by Pic O₂, or they can be user-replaceable, exemplified by the OxySpy model.





SALE AND DELIVERY

Your safety and satisfaction are our top priority

At Nrc International, our employees are dedicated to providing excellent customer service. We pride ourselves on going the extra mile to assist our customers and meet their needs. If you have any questions or concerns, please do not hesitate to reach out to us we are always available to help.



Customer Service

Tel: +49(0)2261/500414
Fax: +49(0)2261/500415
E-mail: info@nrc-international.com



Direct purchase and delivery information

Scan the QR code or visit :
www.NRC-International.com

Notes



NRC International GmbH
Karlstrasse 21a
51702 Bergneustadt
Germany

Info@nrc-international.com
+49 2261 500414

WWW.NRC-INTERNATIONAL.COM