

RENAL CARE

# NIPROSET™

BLOOD TUBING SET FOR HEMODIALYSIS




 **NIPRO**

# NIPROSET™

## BLOOD TUBING SET FOR HEMODIALYSIS

- Vast range of tubing sets for the most prevalent dialysis machine models
- For use in HD and HDF in either Single Needle or Double Needle
- Different sterilization methods available: Gammaphot, E-beam, or Ethylene Oxide (EtO)
- Not made with DEHP\*



## BLOOD TUBING SET PORTFOLIO

Dialysis machine	Article code	Therapy type	Sets/box	Priming volume (mL)
Surdial 55Plus	A146-V715	HD DN	24	129
	A317-V815	HD DN	24	110
B. Braun Dialog, Dialog Advanced, Dialog+	A165-V733	HD DN	24	164
FMC 4008	A183-V751	HD DN	24	157
FMC 2008/4008	A188-V756	HD DN	22	159
FMC 2008/4008	A289-V800	HD DN	22	163
Gambro AK 95-96-98	A164-V732	HD DN	24	127
Gambro AK 95-100-200	A252-V777	HD DN	24	149
Universal	A114-V684	HD DN	24	155
	A134-V721	HD DN	24	155
	A176-V684	HD DN	24	150
	A186-V754	HD DN	24	145


Sterilization method: EtO

Dialysis machine	Article code	Therapy type	Sets/box	Priming volume (mL)
Surdial X	A363R-V849R*	HD DN	24	129
	A372R-V858R*	HD DN	24	129
	A364R-V850R*	HDF DN	20	132
	A365R-V851R*	HD SN	20	231
	A366R-V852R*	HDF SN	18	231
	A367R-V853R*	HD DN	22	131
	A368R-V854R*	HDF DN	20	133
	A369R-V855R*	HD SN	18	231
	A370R-V856R*	HDF SN	16	231
	A438R-V918R*	HDF DN	20	134
	A446R-V926R*	HDF DN	20	133
	A447R-V927R*	HDF SN	16	231
	A146R-V715R	HD DN	24	126
Surdial 55Plus	A317R-V815R	HD DN	24	107
	A165R-V733R	HD DN	24	155
	A290R-V801R	HD DN	24	151
	A320R-V818R	HD DN	24	144
	A434R-V914R	HDF DN	18	163
	A425R-V907R	HD DN	24	154
	A183R-V751R	HD DN	24	169
	A188R-V756R	HD DN	22	159
	A280R-V791R	HD DN BVM	22	175
	A287R-V798R	HD SN	22	248
	A289R-V800R	HD DN	22	156
	A164R-V732R	HD DN	24	128
	A205TR-V773TR*	HD SN BVS	12	420
B. Braun Dialog, Dialog Advanced, Dialog+	A252R-V777R	HD DN	24	145
	A295R-V806R	HD SN	14	415
	A347TR-V719TR*	HD DN BVS	22	167
	A360TR-V846TR*	HD DN BVS	20	146
	A196R-V764R	HD DN	20	163
	A197R-V765R	HDF DN	18	190
	A414R-V897R	HD SN	16	236
	A134R-V721R	HD DN	24	160
	A176R-V684R	HD DN	24	145
	A179R-V747R	HD DN	24	150
	A186R-V754R	HD DN	24	147
	A316R-V814R	HD DN	24	156
	A329R-V827R	HD DN	24	126
	A427R-V909R	HD DN	22	156
Universal				

Sterilization method: Gamma ray or E-beam



Drip chamber



Pillow




Color-coded clamps




Transducer protector



Recirculator



Medication port



Infusion set

Not all products are registered and approved for sale in all countries or regions. Please contact your local Nipro representative for product availability and information.

\*Not made with DEHP



Nipro Renal Care is part of Nipro Corporation Japan, a leading global healthcare company established in 1954. With over 29.000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro Renal Care is a global market leader with over 5 decades providing renal solutions for dialysis and dialysis-related treatment. We specialize in developing dialysis machines, water treatment systems, and a comprehensive portfolio of disposable medical equipment.


In order to address the needs of patients, healthcare professionals, and procurement managers alike, Nipro Renal Care is driven by innovation and patient safety to offer the highest quality products that optimize time, effort, and costs.

## BECAUSE EVERY LIFE DESERVES AFFORDABLE CARE



[www.nipro-group.com/en-en/our-company/our-locations](http://www.nipro-group.com/en-en/our-company/our-locations)

Please contact your local representative for more information.



**Nipro Australia Pty Ltd** : Suite 2.02, Level 2, 657 Pacific Hwy, St Leonards NSW 2065 | T: +61 2 9437 9000 or 1800 451 737 | [info@niproaustralia.com.au](mailto:info@niproaustralia.com.au) | [niproaustralia.com.au](http://niproaustralia.com.au)

**Nipro Medical Europe** : European Headquarters, Blokhuisstraat 42, 2800 Mechelen, Belgium | T: +32 (0)15 263 500 | [www.nipro-group.com](http://www.nipro-group.com)

**Nipro Corporation Japan** : 3-9-3 Honjo-nishi, Kita-ku, Osaka 531-8510, JAPAN | T: +81 6 6372 2331 | [www.nipro.co.jp](http://www.nipro.co.jp)