



SECTION A-A'

NO	DESCRIPTION	REASON	DATE	CHECK
⚠				

REMARK

- NOTE
1. LEAKAGE TEST METHOD : VACUUM SUCKING TEST AT -60cmHg FOR 1MIN/WATER
 2. DOSAGE : 0.20ml±0.02ml / STROKE
 3. FILL CAPACITY : 101.00ml±1.0ml(WATER)
 4. O.F.C. : 106.05ml±1.0ml(WATER)
 5. WEIGHT : 100.95g (PART 00.0±5%, ASS'Y 00.0±10%)
 6. General terms and conditions of quality to comply with Yonwoo standard

9	BOTTOM	PP	EXISTING MOLD	BT195
8	PISTON	PE	EXISTING MOLD	PT192
7	OUTER BOTTLE	PETG	EXISTING MOLD	BOB92
6	INNER BOTTLE	PP	EXISTING MOLD	BOB91
5	PUMP ASS'Y	-	EXISTING MOLD	-
4	SCREW CAP	PP	EXISTING MOLD	SC305
3	SHOULDER	ABS	EXISTING MOLD	PL496
2	ACTUATOR	PP	EXISTING MOLD	BU343
1	OVER CAP	PETG	EXISTING MOLD	OCD17
NO	ITEM	MAT'L	REMARK	PART NO.

TITLE. AS0100-DUL-CC FULL ASS'Y

DWG NO. Q2614-0 SCALE 1 : 1

Drawn by Jung, Woo Jae Checked by [Redacted] Approved by [Redacted] DATE 2021.09.07

REF. TABLE [Redacted]

Engineering 2021.09.07 BongHun

DECIMALS UNIT: mm

Unspecified Tolerance
 0 < 5 ±0.05
 5 < 20 ±0.10
 20 < 40 ±0.20
 40 < 80 ±0.30
 80 < 160 ±0.40
 160 < 320 ±0.50
 Angular ±0.5°

instockpack.