



SECTION A-A'

NO	DESCRIPTION	REASON	DATE	CHECK
△				

REMARK
NOTE
1. LEAKAGE TEST METHOD : VACUUM SUCKING TEST AT -60cmHg FOR 1MIN/WATER
2. DOSAGE: 0.20ml ± 0.02ml / STROKE
3. FILL CAPACITY: 16.00ml ± 1.0ml(WATER)
4. O.F.C.: 18.50ml ± 1.0ml(WATER)
5. WEIGHT : 24.00g (PART 00.0±5%, ASS'Y 00.0±10%)
6. General terms and conditions of quality to comply with Yonwoo standard

NO	ITEM	MAT'L	REMARK	PART NO.
8	BOTTOM	PP	EXISTING MOLD	BT028
7	PISTON	PE	EXISTING MOLD	PT029
6	BOTTLE	PP	EXISTING MOLD	BO094
5	PUMP ASS'Y	-	EXISTING MOLD	-
4	SNAP ON CAP	PP	EXISTING MOLD	SO008
3	PUMP CAP	PP	EXISTING MOLD	PL032
2	ACTUATOR	PP	EXISTING MOLD	BU102
1	OVER CAP	PP	EXISTING MOLD	OC096

TITLE.	AS0015-S/PP FULL ASS'Y			
DWG NO.	Q1898-0	SCALE	1 : 1	
Drawn by	Checked by	Checked by	Approved by	
Nam, Mi Hye				Unspecified Tolerance 0 < 5 - ±0.05 5 < 20 - ±0.10 20 < 40 - ±0.20 40 < 80 - ±0.30 80 < 160 - ±0.40 160 < 320 - ±0.50 Angular - ±0.5°
DATE	2017.12.07	DECIMALS UNIT: mm		