



SECTION A-A'

NO	DESCRIPTION	REVISION	DATE	CHECK
△				

REMARK

NOTE

1. LEAKAGE TEST METHOD : VACUUM SUCKING TEST AT -60cmHg FOR 1MIN/WATER
2. DOSAGE: 0.20ml± 0.02ml / STROKE
3. FILL CAPACITY: 51.0ml±1.0ml(WATER)
4. O.F.C.: 62.5ml±1.0ml(WATER)
5. WEIGHT : 41.00g (PART 00.0±5%, ASS'Y 00.0±10%)
6. General terms and conditions of quality to comply with Yonwoo standard

8	PISTON	PE	EXISTING MOLD	PT147
7	BOTTLE	PP	EXISTING MOLD	BO797
6	HOUSING	PP	EXISTING MOLD	HU156
5	SHOULDER CAP	PP	EXISTING MOLD	PL333
4	PUMP ASS'Y	-	EXISTING MOLD	-
3	COLLAR	PP	EXISTING MOLD	PL287
2	ACTUATOR	PP	EXISTING MOLD	BU323
1	OVER CAP	PETG	EXISTING MOLD	OC294
NO	ITEM	MAT'L	REMARK	PART NO.

TITLE. AS0050-PN FULL ASS'Y

DWG NO. Q1915-0 SCALE 1 : 1

Drawn by Nam, Mi Hye

Checked by

Checked by

Approved by **Engineering** 2017.12.22 BongHun

MAT'L REF.TABLE DATE 2017.12.22

DECIMALS UNIT: mm

Unspecified Tolerance  
 0 < 5 - ±0.05  
 5 < 20 - ±0.10  
 20 < 40 - ±0.20  
 40 < 80 - ±0.30  
 80 < 160 - ±0.40  
 160 < 320 - ±0.50  
 Angular - ±0.5°

**instockpack.**