


ComLab 3-Way Bench Mounted Laboratory Mixing Set with 2 Bib Outlets and 1 Swivel/ Fixed Gooseneck Outlet - Ceramic Disc - Water

Enware Commercial Laboratory outlets are supplied with a tapered tube nozzle and are assembled for each specific media type. The durable powder coat finish gives significant protection against corrosion, discolouration and surface damage. Slimline handles are ergonomically designed for precise feel and fingertip control. Handles are colour coded for easy service identification and use complimentary buttons.



KEY FEATURES

- 3/4-turn ceramic disc headworks provides precise flow control
- Durable white epoxy powder coat finish with coloured metal handles
- 2 x Fixed bib outlets (hot and cold)
- Swivel high-rise gooseneck spout with integral fixing option
- 1/2" BSP male inlet tails with backnuts for bench mounting
- Tube nozzle outlet with flow control and dual check valve. 3/8" BSP FI when removed
- Aerator option
- Green handle and indicator button identification to DIN EN13792
- Indicator buttons for compatible media included
- Compatible media: Cold Water, Hot Water, Chilled Water
- Lead free brass

PRODUCT CODE

LCN108

ComLab 3-Way Bench Mounted Laboratory Mixing Set with 2 Bib Outlets and 1 Swivel/Fixed Gooseneck Outlet - Ceramic Disc - Water



Due to ongoing Research and Development, specifications may change without notice.

Component specifications may change on some export models.

Refer to warranty statement for warranty details - www.enware.com.au/warranty.

Products are to be installed in accordance with the Plumbing Code of Australia and AS/NZS3500.

Reference should also be made to the Australasian Health Facility Guidelines (AHFG), ABCB and Local Government regulations when considering the choice of, and the installation of these products.

ComLab 3-Way Bench Mounted Laboratory Mixing Set with 2 Bib Outlets and 1 Swivel/ Fixed Gooseneck Outlet - Ceramic Disc - Water



TECHNICAL DATA

WELS Category	Tapware
Standard WELS	5 Stars
Standard WELS Flow Rate	6 L/min
Finish	Powder Coating
Valve Type	Ceramic Disc
Water Options	Cold, Hot, Chilled
Maximum Working Temperature	90°C
Maximum Recommended Pressure	700kPa
Maximum Static Pressure	1000kPa
Primary Flow Rate	6 L/min
Inlet Connection	1/2" BSP Male
Outlet Connection	Tube Nozzle, 3/8"BSP Female (without nozzle)

Tube Nozzle may cause splashing and is designed to have tubing fitted when in use.

Aerator is available as an option.

Water used in the laboratory may not be potable, additional signage may be required.



INSTALLATION DETAIL
(BENCH TOP CUT OUT)

OPTIONAL ACCESSORIES

LCN172-5	Aerator Kit - 5.7 L/min Laminar Flow to suit LCN Products
LCN175	Spare Parts: Laboratory - LCN Adaptor - 1/2" BSP M Outlet

SPARE PARTS

LCN183	ComLab: Stuffing Box Assembly - 3/4 Turn Ceramic Disc - LCN Products
LCN182	ComLab: Button Replacement - Water - LCN Products
LCN187	Dual Check Valve - LCN Nozzle
LCN190	Flow Control - LCN Tube Nozzle – 8 L/min
LCN191-GREEN	ComLab: Handle (Green)
LCN171	Spare Parts: ComLab - Tube Nozzle - Dual Check Valve - LCN Water Products
LCN193	ComLab: Indicator Mount
FS734	Spare Parts: Washer Service Kit plus O-Ring
LCN186	ComLab: Base Flange for Wall/Bench