

ESSENTIAL WORKPLACE RESPONSE KIT - REFILL

PART NUMBER: FAEWR

CARDBOARD REFILL



DESCRIPTION:

- National workplace first aid kit – Low risk.
- Designed in accordance with Safe Work Australia’s First Aid Code of Practice (July 2019).
- Packaged in a white cardboard box, ideal for refilling a complete kit or replenish a first aid kit as required.
- This a complete refill pack to suit:
 - FAEWM Kit
 - FAEWT Kit
 - FAEWS Kit

LOW RISK KIT OVERVIEW:

Count: 1-25 Person Kit	Sold As: Each	Pieces: 309 (Contents Below)

CONTENTS & QUANTITY

Accident Report Notebook With Pencil	1	Non-Adherent Dressing Sterile, 10 X 10cm	2
Adhesive Strips Plastic, 72 X 19mm	50	Non-Adherent Dressing Sterile, 10 X 7.5cm	3
Alcohol Wipe Sachet	4	Non-Adherent Dressing Sterile, 5 X 5cm	6
Antiseptic Liquid/Sting Relief Spray Bottle, 50ml	1	Paper Tape Hypoallergenic White, 2.5cm X 5m	1
Bowl Galipot Clear, 60ml	1	Plastic Bag Resealable Large	1
Cold Pack Instant Large, 160g	1	Plastic Bag Resealable Medium	1
Combine Dressing Sterile, 20 X 20cm	1	Plastic Bag Resealable Small	1
Conforming Bandage White, 5cm	3	Povidone Iodine Swabs	6
Conforming Bandage White, 7.5cm	3	Resuscitation Face Shield Disposable With Non-Return Valve	1
Cotton Tips Double Ended, 7.5cm	100	Safety Pins, Assorted	12
CPR Flowchart Wallet Size	1	Scissors Stainless Steel, 12.5cm	1
Crepe Bandage Heavy Weight Brown, 10cm	1	Sharps Container Plastic Clear, 200 MI	1
Disposable Gloves Large, (5 Pairs)	10	Splinter Forceps Stainless Steel, 12.5cm	1
Emergency First Aid Information Booklet Basic DL Size	1	Splinter Probe Disposable, 5pcs/Sheet	10
Emergency Shock Blanket Silver, 160 X 210cm	1	Triangular Bandage Cotton Large White	2
Eye Pad Sterile	4	Wound Cleansing Wipe Non-Sting Sachet	10
Eye Wash Solution, 20ml Ampoule	8	Wound Closure Strips, 3 X 75mm, 5pcs/Sheet	25
Gauze Swab, 7.5 X 7.5cm, Pack Of 5	25	Wound Dressing No.13 Sterile	1
Hydrogel Burns Gel, 3.5g Sachet	5	Wound Dressing No.14 Sterile	1
Kidney Tray Plastic Clear, 500ml	1	Wound Dressing No.15 Sterile	1

*All contents maybe subject to change