

 **VAPOTHERM®**

Contract #: PP-NS-1730
Category: Active Humidification Devices and Accessories
Start date: November 1, 2023
End date: April 30, 2024
PO deadline: June 21, 2024
Product: HVT 2.0 with Transfer Stand

- Members **save up to 25% off best national tier!**
- No minimum order required. Maximum savings for 10+ units.
- Pricing includes HVT 2.0 Unit with Transfer Stand.
- Lock in your discounts on **stockd.com** by April 30, 2024. Orders are confirmed and fulfilled directly with Vapotherm; Purchase Orders (POs) must be issued by June 21, 2024.
- Exclusive offer for Premier GPO members only.*

Proven Therapy. Simply Delivered.




**CLICK HERE TO LOCK IN
YOUR DISCOUNTS TODAY!**

** Members from digital sales restricted classes of trade (COTs) - Affiliate Physician(s) and Oncology Centers - must contact Vapotherm directly to participate; all other members participate by locking in discounts on stockd®.*

learn more at **stockd.com**

Why HVT 2.0?

HVT 2.0 is simple to use and fast to set up.
Built to work for you and your patients.



Vapotherm HVT 2.0 uses our revolutionary high velocity therapy and is **ready for use in approximately 5 minutes**. It requires minimal patient coaching, and staff intervention.



A **large, intuitive screen** gives you clinically relevant feedback from a distance, which means you can make informed decisions quickly, and with confidence.

HVT 2.0 maintains patient safety and infection control with a disposable patient circuit that is **simple to change and set-up between patients**. It provides on-screen troubleshooting support, and comprehensive alarms that facilitate patient safety and ease of care.

Untethered Respiratory Support

Seamless mobility and connectivity for your patients.
Anytime. Anywhere.



Our solution is **powered by AC with an internal battery** to add safety in the event of a power failure, a **transfer battery** for additional power needs, and has an **internal blower** that eliminates the need for air tanks during patient transfer.

Mounted on our ergonomically designed roll stand and engineered to hold two O₂ cylinders, HVT 2.0 allows you to safely maintain therapy, when and where you need it throughout the hospital - untethered from the wall with optimal mobility.

learn more at stockd.com