

Molift Partner 255

The Molift Partner 255 is a perfect solution for hospitals and care facilities. With a high lifting capacity, it can lift from the floor up. The Molift Partner 255 has a vertical lifting pattern which keeps the user at an equidistant length from the column during the entire lift, creating a safe and stable position.



Properties

- Secure and comfortable
- Easy to maneuver
- Lifts from the floor
- Lifts users vertically
- Environmentally friendly batteries

Molift Partner 255

When the lifter is used with the Molift 4-point suspension and Molift slings, this provides comfortable and correct positioning for lying to sitting and vice versa. Low weight and superior wheel positioning makes the Molift Partner 255 incredibly easy to manoeuvre.

The lifter combined with the low chassis and electric leg positioning means that the Molift Partner 255 can be used in most situations. The Molift Partner 255 is equipped with a built-in automatic service notification. It has electronic/manual emergency lowering functions and emergency stop. Molift Partner 255 is available with an extra low chassis, only 6 cm.



All passive 4 point Molift slings can be used with the Molift Partner 255.

Technical facts

Item no.

xxx

Lifting capacity

255 kg

Total weight

44 kg

Lifting height

min 400 mm -
max 1680 mm

Chassis width

700 mm

Total Width

68.5 cm inner/124.46 cm outer

Leg height

115 mm / 60 mm low chassis

Battery

NiMH 26.4 V - 2.2 Ah

Accessories

**Molift Partner 255 Scale**

The scale complies with all applicable technical standard and is approved class III scale. The scale can calculate Body Mass Index.

Item no: 1331002

**Molift Stretcher**

A scoop carrier specially adapted to users with back injuries. It is X-ray compatible and can be suspended directly from the standard 4-point suspension.

Item no: 2150103

**Molift Partner Walking Arms**

Molift Partner 255 is available with adjustable support arms that can be used with the Molift ambulating vest for walking/training. The arms fold along the column when not in use.

Item No.: 1311052