Molift Mover 300

With its high lifting capacity Molift Mover 300 is the perfect solution for hospitals and care facilities. It has tall lifting height, and can therefore lift from the floor up. The lowest weight in its category in the market combined with the superior wheel positioning makes the Molift Mover 300 incredibly easy to maneuver. The low chassis, only 11 cm, together with the electric leg positioning provides for easy access.





Molift Mover 300

The Molift Mover 300 is easy to clean and disinfect. Molift Mover 300 is supplied with 4-point suspension as standard. Used together with a Molift Sling, it provides comfortable and correct positioning from lying- to- sitting and vice versa. Molift Mover 300 can also be supplied with alternative suspensions to produce a quick coupling system that facilitates quick and simple switches.

Molift Mover 300 is equipped with automatic service notification. This consists of integrated software which calculates when the lifter requires a service, as well as produces statistics on the user's lifting patterns. It is also equipped with electronic and manual emergency lowering and emergency stop.



Molift Mover 300 can be used with all passive Molift slings designed for 4-point suspension.

Technical facts

Item no.

15000

LIfting capacity 300 kg

Total weight 39.3 kg

Lifting height min 400 mm max 1750 mm Chassis width 700 mm

Leg height 115 mm

Battery

NiMH 26.4 V - 2.2 Ah

Accessories



Molift Mover 300 Scales
The scale complies with all
applicable technical standard
and is approved class III scale.
The scale can calculate Body
Mass Index.

Item no: 1431012



Molift Stretcher
A scoop carrier specially adapted to users with back injuries. It is X-Ray, CT Scan, MRI compatible and can be suspended directly from the standard 4-point suspension.

Item no: 2150103 (standard) og 2150108 (MRI compatible)



Molift Mover Walking Arms
Molift Mover 300 is available with
adjustable support arms that can
be used with the Molift ambulating vest for walking/training. The
arms fold along the column when
not in use.
Item No.: 1411002

LIFTABILITY molift