

K199 THE SPOONER STRAP

Patella Stabilizing Strapping System

Recommended LCodes:

Custom - L2800

Non-Custom - L1820

The K199, also known as the "Spoooner," is a product that can help with patella-femoral tracking issues while wearing a rigid knee brace after a knee ligament injury or surgery. This product can be used with any brand of rigid brace and can be strapped to the brace or worn directly on the skin to alleviate patella soreness caused by mal-tracking.

When patella subluxation occurs in conjunction with multiple instabilities of the knee, the K199 strapping system can be fitted to most braces on the market today to simultaneously control patella instability problems and migration issues, providing a comprehensive solution for knee support and stability.

FEATURES:

- Easily allows for retrofitting to most rigid braces on the market today without modifying existing products.
- "The Spooning" effect captures and controls the patella.
- Cost-effective way to prevent functional brace migration.
- The Spooner is the practical solution to common treatment issues with rigid braces

INDICATIONS:

- Controls medial patella Meniscus Ligament (MPML) and Medial Femoral Ligament (MPFL)
- Assures proper alignment while wearing most rigid functional braces.

Washing Instructions:

Hand wash in cold water with mild detergent.
NO chlorine bleach. Air dry. Do NOT put in dryer.



1850 Diplomat Drive, Suite #100
Farmers Branch, TX 75234
PH#: 214-638-6422 FX#: 214-638-6425

