

KC64CT-NOS KneeMate™

COLD THERAPY BRACE

Recommended LCodes:
L1832/L1833

FEATURES:

- Continuous cold therapy without removing the brace
- Multi-use cold pads (not included) attach to brace interior
- Hinges can be inverted for future wear as rehab brace
- Sharkskin anti-migration strips
- Fits universal left or right
- Upper & lower circumferential straps with elastic segments
- Breathable, open-cell KoolFlex material
- Hinges with Flexion and Extension stops

INDICATIONS:

- Post-operative
- Arthroscopic Procedures
- Post-Trauma
- Chronic Pain
- Physical Therapy
- Medial/Lateral instabilities where cryotherapy is indicated



Washing Instructions:

Hand wash in cold water with mild detergent.
NO chlorine bleach. Air dry. Do NOT put in dryer.



SIZE	KNEE CIRCUMFERENCE	6" ABOVE KNEE
X-Small	12" - 13"	16" - 19"
Small	13" - 14"	17" - 20"
Medium	14" - 15"	18" - 21"
Large	15" - 16"	19" - 22"
X-Large	16" - 18"	20" - 24"
2X-Large	18" - 20"	22" - 26"
3X-Large	20" - 22"	24" - 28"
4X-Large	22" - 24"	26" - 30"
Custom	24" & UP - PLEASE PROVIDE MEASUREMENTS	

1850 Diplomat Drive, Suite #100

Farmers Branch, TX 75234

PH#: 214-638-6422 FX#: 214-638-6425

