

KneeMate™ Series

PDAC Approved: KC64-NOS & KC64-MP - L1832/L1833
 PDAC Approved: KC64-PC - L1820
 Recommended LCode: KC64-SL - L1810

FEATURES:

- Easy anterior hook and loop closure with 5 inches of adjustable hook & loop
- Adjustable hook & loop closures with easy open pull tabs
- Comfortable large popliteal opening
- Hinges enclosed in neoprene pocket or open
- Added elastic strap segment for improved fit and comfort
- Special cordura stress covers over center proximal and distal hinges

Product Series:

K64-PC	With steel, polycentric hinges
K64-MP	With multi-position hinges
K64-HT	With high tech, polycentric hinges
K64-E	With spiral stays with dipped caps
K64-NOS	With range of motion hinge

INDICATIONS:

- Medial & Lateral collateral instability
- Arthritis
- Post-Op management
- Mild to moderate anterior and posterior cruciate instability



Washing Instructions:

Hand wash in cold water with mild detergent.
 NO chlorine bleach. Air dry. Do NOT put in dryer.



KC64-NOS



KC64-MP



KC64-SL



KC64-PC

KoolFlex™ versions:

KC64-E, KC64-SL, KC64-PC, KC64-MP, KC64-HT, KC64-NOS

SIZE	KNEE CIRCUMFERENCE	6" ABOVE KNEE
X-Small	12" - 13"	16" - 19"
Small	13" - 14"	17" - 20"
Medium	14" - 15"	18" - 21"
Large	15" - 16"	19" - 22"
X-Large	16" - 18"	20" - 24"
2X-Large	18" - 20"	22" - 26"
3X-Large	20" - 22"	24" - 28"
4X-Large	22" - 24"	26" - 30"
Custom	24" & UP - PLEASE PROVIDE MEASUREMENTS	

1850 Diplomat Drive, Suite #100

Farmers Branch, TX 75234

PH#: 214-638-6422 FX#: 214-638-6425

