

# **OAB: PULL-ON OA BRACE**

Recommended LCodes: L1843/1851; L1844 for Customs

### **Introducing the OA Brace**

Exceptional comfort and simplicity. Precise support for the osteoarthritis knee patient.

Innovative construction helps patients comply with your treatment plan and enjoy the freedom of a healthy lifestyle.

## **FEATURES:**

- Compressive buttressed knee sleeve, assuring correct joint alignment
- Soft construction for stability without "over bracing"
- Simple adjustment and fit with proximal and distal closures
- Strut paddles for rotation control and load
- Pre-set "Ezy Hinge" for optimal counter force and extension control
- Low profile, lightweight design for discreet support beneath clothing
- · Quality, American-made craftsmanship

## **INDICATIONS:**

- Mild to moderate levels of unicompartmental osteoarthritis (OA) of the knee.
- · Post meniscus repairs
- · Micro-fracture surgery



#### Washing Instructions:

Hand wash in cold water with mild detergent. NO chlorine bleach. Air dry. Do NOT put in dryer.







#### **SIZING**

Measure knee circumference and 6" above the center of the knee with leg extended

| SIZE     | KNEE<br>CIRCUMFERENCE | 6" ABOVE KNEE |
|----------|-----------------------|---------------|
| X-Small  | 12" - 13"             | 16" - 19"     |
| Small    | 13" - 14"             | 17" - 20"     |
| Medium   | 14" - 15"             | 18" - 21"     |
| Large    | 15" - 16"             | 19" - 22"     |
| X-Large  | 16" - 18"             | 20" - 24"     |
| 2X-Large | 18" - 20"             | 22" - 26"     |
| 3X-Large | 20" - 22"             | 24" - 28"     |
| 4X-Large | 22" - 24"             | 26" - 30"     |
| Custom   | 24" + Up              | +30"          |

1850 Diplomat Drive, Suite #100 Farmers Branch, TX 75234 PH#: 214-638-6422 FX#: 214-638-6425





# OAB: PULL-ON OA BRACE

Suggested LCodes: L1843/1851; L1844 for Customs

## **ORDERING:**

When placing an order, select either OABL or OABR. OABL products are designed for left knee medial OA or right knee lateral OA patients, with the hinge on the left, while OABR products are intended for right knee medial OA or left knee lateral OA patients, with the hinge on the right. Both options are available in X-Small through 4X-Large sizes.

Example: A Large Right Medial/Left Lateral OA Brace part number is: OABR-LG.

| ITEM CODE | DESCRIPTION   |
|-----------|---|
| OABL-XS   | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, XS  |
| OABL-S    | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, S   |
| OABL-M    | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, M   |
| OABL-LG   | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, LG  |
| OABL-XL   | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, XL  |
| OABL-2XL  | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, 2XL |
| OABL-3XL  | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, 3XL |
| OABL-4XL  | OA Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, Blue, 4XL |

| ITEM CODE | DESCRIPTION   |
|-----------|---|
| OABR-XS   | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, XS  |
| OABR-S    | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, S   |
| OABR-M    | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, M   |
| OABR-LG   | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, LG  |
| OABR-XL   | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, XL  |
| OABR-2XL  | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, 2XL |
| OABR-3XL  | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, 3XL |
| OABR-4XL  | OA Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, Blue, 4XL |



1850 Diplomat Drive, Suite #100 Farmers Branch, TX 75234 PH#: 214-638-6422 FX#: 214-638-6425

