

FLEX® OTS

Flex is a viable post-op/functional rehab solution offering more versatility in both shape and function than most other options. Its heat-moldable, carbon composite frame accommodates edema and/or changes in leg shape during recovery from surgery, while its Accutrac® range-of-motion hinges provide flexion and extension control. Gel Fit condylar pads improve suspension and are comfortable for tender knees recovering from trauma.







INDICATIONS

- · ACL, MCL, LCL, PCL, rotary and combined instabilities
- Long-term use for post-op, rehab and functional support
- Recommended for low to medium contact/impact
- Sport model for lean, athletic legs

PART#	LEG	SIZE	MEASURING ¹
FLEX-SR	Right	Small	Caliper: 3.5-4"
FLEX-SL	Left	Small	Circumference: 11-12.625"
FLEX-MR	Right	Medium	Caliper: 4-4.625"
FLEX-ML	Left	Medium	Circumference: 12.625-14.5"
FLEX-LR	Right	Large	Caliper: 4.625-5.125"
FLEX-LL	Left	Large	Circumference: 14.5-16"
FLEX-XLR	Right	XLarge	Caliper: 5.125-5.625"
FLEX-XLL	Left	XLarge	Circumference: 16-17.625"

¹ Caliper: M/L at joint width in a standard, weight-bearing position Circumference: just below the distal border of the patella

Flex Sport

PART#	LEG	SIZE	MEASURING ¹
FLEXSP-SR	Right	Small	Caliper: 3.5-4"
FLEXSP-SL	Left	Small	Circumference: 11-12.625"
FLEXSP-MR	Right	Medium	Caliper: 4-4.625"
FLEXSP-ML	Left	Medium	Circumference: 12.625-14.5"
FLEXSP-LR	Right	Large	Caliper: 4.625-5.125"
FLEXSP-LL	Left	Large	Circumference: 14.5-16"
FLEXSP-XLR	Right	XLarge	Caliper: 5.125-5.625"
FLEXSP-XLL	Left	XLarge	Circumference: 16-17.625"

¹ Caliper: M/L at joint width in a standard, weight-bearing position Circumference: just below the distal border of the patella