



OBSERVA SERIES

PRIMALED

The first and only examination lamp with the possibility to choose an ambient light in addition to traditional concentrated light. The Ergo-Spring balancing system makes PRIMALED very easy to handle and stable. PRIMALED is ideal for any type of installation, from the outpatient department to the intensive care unit.

Battery Group



Fixing Systems





PRIMA-FIX
Joint Arm



PRIMA-FLEX
Gooseneck Arm

Performances

PRIMALED

| | |
|--|----------------|
| Light intensity at 0,50 m distance (Ec) | 105.000 lux |
| Light head diameter | 19,5 cm |
| Color temperature (2 selections) | 4.000 - 4.500K |
| Color rendering Index (CRI) | 95 Ra |
| d10 light field diameter where illuminance reached 10% of Ec | 150 mm |
| Average LED life | 60.000 hours |
| Control of the illuminance | 20 - 100 % |
| Electrical absorption | 10,5 W - 20 VA |

MAIN FEATURES



ACCESSORIES



Regulation (EU) 2017/745 of the European Parliament and of the Council, of 05th April 2017 – Norm IEC 60601-2-41

Contact **THE LAMPHOUSE MEDICAL**
info@thelamphouse.co.za
011 608 2565