





Designed by Edie Eckman

# What you will need:

RED HEART® Soft®: 1 ball 9388 Wheat

Susan Bates® Knitting Needles: 5 mm [US 8]

Cable needle, yarn needle

GAUGE: 19 sts = 4" (10 cm); 23 rows = 4" (10 cm) in Stockinette st (knit on right side, purl on wrong side). CHECK YOUR GAUGE. Use any size needles to obtain the gauge.



RED HEART® Soft®, Art. E728 available in 5 oz (141 g), 256 yd (234 m) balls



# Basic Cables Square

For Knit Your Cables Afghan

Use this square for making the one-color Knit Your Cables Afghan. Follow our blog by searching RHK0502-019189M for hints on knitting each square of this afghan. Then, if you feel creative, use the stitch pattern for your own unique knit project. Square measures 10" (25.5 cm) square.

#### **Special Abbreviation**

**2/2 LC (2 over 2 Left Cross):** Slip next 2 sts to cable needle and hold in front, k2, k2 from cable needle.

### **Special Technique**

Cable Cast On = Place a slip knot on left needle, insert right needle into stitch on left needle as if to knit, pull loop through and place on left needle (2 sts on left needle), \*insert needle between first 2 stitches on left needle, pull loop through and place on left needle; repeat from \* until the required number of stitches are on the left needle.

## **SOUARE B**

Cable cast on 50 sts.

**Rows 1 and 5 (right side):** P7, [k4, p4] 5 times, p3.

**Rows 2 and 4:** K3, [k4, p4] 5 times, k7. **Row 3:** P7, [2/2 LC, p4] 5 times, p3.

Row 6: Repeat Row 2.

Repeat Rows 1-6 ten more times. Bind off in pattern on right side.


## **Finishing**

Block piece to approximately 10" (25.5 cm) square.

Weave in ends.

#### **ABBREVIATIONS**

 $\mathbf{k} = \text{knit}; \ \mathbf{p} = \text{purl}; \ \mathbf{st(s)} = \text{stitch(es)}; \ [\ ] = \text{work}$  directions in brackets the number of times specified.



8-st repeat

**KEY** 

- □ k on RS, p on WS
- □ p on RS, k on WS
- ☑ k tbl

- △ 2/1 LPC
- **Z** 2/1 RPC
- ≥ 1/2 LPC
- 1/2 RPC
- 2/2 RC
- ≥ 2/2 LC
- □ 3/1 LPC
- **⊘** 3/1 RPC
- 3/3 RC 3/3 LC

**SHOP KIT**