

## dNTP Mix (10 mM each)

PRODUCT	SIZE	SKU
dNTP Mix (10 mM each)	0.5 ml	TRAXSKU1019
	1 ml	TRAXSKU1020
	2.5 ml	TRAXSKU1021

### STORAGE & HANDLING

Store at -20°C. Thaw and mix properly before use. Avoid excessive freeze and thawing cycles, prepare aliquots if necessary.

### DESCRIPTION

dNTPs are an essential component of every PCR reaction. dATP, dTTP, dCTP and dGTP are supplied in a convenient ready-to-solution with 10 mM of each nucleotide that can be used directly for most PCR protocols without further dilution.

### FEATURES

- Ultrapure: ≥99% as determined by HPLC
- Supplied in ultrapure dH<sub>2</sub>O
- Nuclease-free
- Highly stable
- Compatible with all Procomcure PCR reagents e.g. VitaTaq® DNA Polymerase