


Analysis of polychlorinated dibenzodioxins and furans


BPX5

Column part number	054114		
Phase	BPX5	Rate 3	10°C/min to 320°C
Column	50 m x 0.22 mm x 0.25 µm	Final temperature	320°C
Initial temperature	80°C	Detector	MS
Rate 1	10°C/min to 240°C	Carrier gas	He
Rate 2	2°C/min to 280°C		

Single ion chromatograms of TCDF isomers



Single ion chromatograms of TCDD isomers



Acknowledgements: Dr P.Ambridge, Dr A.Fernandes and C.Brook at AEA Technology, Harwell, UK.

For more information visit www.trajanscimed.com or contact techsupport@trajanscimed.com