


Analysis of organochlorine pesticides (USEPA 8081)

BPX50

Column part number		054751	
Phase	BPX50	Final temperature	320 °C
Column	30 m x 0.25 mm x 0.25 µm	Detector	MS
Initial temperature	50°C, 1 min	Injection mode	Split (40:1), 280°C
Rate 1	20°C/min to 170°C	Carrier gas	He, 35 cm/sec
Rate 2	8°C/min to 320°C	Injection volume	1 µL


 For more information visit www.trajanscimed.com or contact techsupport@trajanscimed.com