


Analysis of alkanes

BP BOND U

Column part number	0571237		
Phase	BP BOND U	Final temperature	250°C, 1 min
Column	25 m x 0.32 mm x 7 µm	Detector	FID, 250°C
Sample	Alkanes 1000 ppm each	Carrier gas	He
Initial temperature	40°C, 1 min	Carrier gas flow	30 mL/min
Rate	15°C/min	Injection mode	Split/splitless


 For more information visit www.trajanscimed.com or contact techsupport@trajanscimed.com