

Not all syringe part numbers are recognized by the Thermo Scientific TriPlus RSH instrument software. When using the SGE syringe part numbers below, enter the corresponding part number for instrument configuration.

SGE syringe part number	Description	Part number for instrument configuration
000480	0.5BR-RTC/RSH-5.7/0.63C 0.5UL SYRINGE	365A0241
000680	1BR-RTC/RSH-5.7/0.63C 1UL SYRINGE	365B0251
001861	5F-RTC/RSH-5.7/0.63C 5UL SYRINGE	365C0231
001863	5F-RTC/RSH-8.5/0.63C 5UL SYRINGE	365C0221
001865	5F-RTC/RSH-5.7/0.47C 5UL SYRINGE	365C0251
001867	5F-RTC/RSH-8.5/0.47C 5UL SYRINGE	365C0241
001871	5R-RTC/RSH-5.7/0.63C 5UL SYRINGE	365C0231
001873	5R-RTC/RSH-8.5/0.63C 5UL SYRINGE	365C0221
001875	5R-RTC/RSH-5.7/0.47C 5UL SYRINGE	365C0251
001877	5R-RTC/RSH-8.5/0.47C 5UL SYRINGE	365C0241
002790	10R-RTC/RSH-GT-MS1 10UL DIAMOND MS SYRINGE	365D0311
002852	10F-RTC/RSH-GT-5.7/0.63(0.24)C 10UL SYRINGE	365D0311
002854	10F-RTC/RSH-GT-8.5/0.72LC 10UL SYRINGE	365D0261
002861	10F-RTC/RSH-5.7/0.63C 10UL SYRINGE	365D0271
002862	10F-RTC/RSH-GT-5.7/0.63C 10UL SYRINGE	365D0311
002863	10F-RTC/RSH-8.5/0.63C 10UL SYRINGE	365D0261
002864	10F-RTC/RSH-GT-8.5/0.63C 10UL SYRINGE	365D0301
002865	10F-RTC/RSH-5.7/0.47C 10UL SYRINGE	365D0291
002866	10F-RTC/RSH-GT-5.7/0.47C 10UL SYRINGE	365D0331
002867	10F-RTC/RSH-8.5/0.47C 10UL SYRINGE	365D0281
002868	10F-RTC/RSH-GT-8.5/0.47C 10UL SYRINGE	365D0321
002869	10F-RTC/RSH-8.5/0.47 BV 10UL SYRINGE	365D0281
002871	10R-RTC/RSH-5.7/0.63C 10UL SYRINGE	365D0271
002872	10R-RTC/RSH-GT-5.7/0.63C 10UL SYRINGE	365D0311
002873	10R-RTC/RSH-8.5/0.63C 10UL SYRINGE	365D0261
002874	10R-RTC/RSH-GT-8.5/0.63C 10UL SYRINGE	365D0301
002875	10R-RTC/RSH-5.7/0.47C 10UL SYRINGE	365D0291
002876	10R-RTC/RSH-GT-5.7/0.47C 10UL SYRINGE	365D0331
002877	10R-RTC/RSH-8.5/0.47C 10UL SYRINGE	365D0281
002878	10R-RTC/RSH-GT-8.5/0.47C 10UL SYRINGE	365D0321
002880	10R-RTC/RSH-GT-5.7/0.72LC 10UL SYRINGE	365D0311
003790	25R-RTC/RSH-GT-MS1 25UL DIAMOND MS SYRINGE	365F2441
003852	25F-RTC/RSH-GT-5.7/0.63(0.24)C 25UL SYRINGE	365F2441
003854	25F-RTC/RSH-GT-8.5/0.72LC 25UL SYRINGE	365F2431
003862	25F-RTC/RSH-GT-5.7/0.63C 25UL SYRINGE	365F2441

SGE syringe part number	Description	Part number for instrument configuration
003864	25F-RTC/RSH-GT-8.5/0.63C 25UL SYRINGE	365F2431
003866	25F-RTC/RSH-GT-5.7/0.47C 25UL SYRINGE	365F2461
003868	25F-RTC/RSH-GT-8.5/0.47C 25UL SYRINGE	365F2451
004790	50R-RTC/RSH-GT-MS1 50UL DIAMOND MS SYRINGE	365G2311
004862	50F-RTC/RSH-GT-5.7/0.63C 50UL SYRINGE	365G2311
004864	50F-RTC/RSH-GT-8.5/0.63C 50UL SYRINGE	365G2321
004866	50F-RTC/RSH-GT-5.7/0.47C 50UL SYRINGE	365G2331
004868	50F-RTC/RSH-GT-8.5/0.47C 50UL SYRINGE	365G2321
005790	100R-RTC/RSH-GT-MS1 100UL DIAMOND MS SYRINGE	365H2141
005852	100F-RTC/RSH-GT-5.7/0.63(0.24)C 100UL SYRINGE	365H2141
005854	100F-RTC/RSH-GT-8.5/0.72LC 100UL SYRINGE	365H2171
005862	100F-RTC/RSH-GT-5.7/0.63(0.11)C 100UL SYRINGE	365H2141
005864	100F-RTC/RSH-GT-8.5/0.63(0.11)C 100UL SYRINGE	365H2171
005866	100F-RTC/RSH-GT-5.7/0.47(0.11)C 100UL SYRINGE	365H2161
005868	100F-RTC/RSH-GT-8.5/0.47(0.11)C 100UL SYRINGE	365H2151
005890	100F-RTC/RSH-GT-5.7/0.63(0.32)H 100UL SYRINGE	365H2181
005891	100F-RTC/RSH-GT-8.5/0.63(0.32)H 100UL SYRINGE	365H2171
006854	250F-RTC/RSH-GT-8.5/0.72LC 250UL SYRINGE	365I2341
006862	250F-RTC/RSH-GT-5.7/0.63C 250UL SYRINGE	365I2351
006864	250F-RTC/RSH-GT-8.5/0.63C 250UL SYRINGE	365I2341
006866	250F-RTC/RSH-GT-5.7/0.47C 250UL SYRINGE	365I2331
006868	250F-RTC/RSH-GT-8.5/0.47C 250UL SYRINGE	365I2321
006890	250F-RTC/RSH-GT-5.7/0.63H 250UL SYRINGE	365I2351
006891	250F-RTC/RSH-GT-8.5/0.63H 250UL SYRINGE	365I2341
007862	500F-RTC/RSH-GT-5.7/0.63C 500UL SYRINGE	365J2441
007864	500F-RTC/RSH-GT-8.5/0.63C 500UL SYRINGE	365J2431
007866	500F-RTC/RSH-GT-5.7/0.47C 500UL SYRINGE	365J2421
007868	500F-RTC/RSH-GT-8.5/0.47C 500UL SYRINGE	365J2411
007890	500F-RTC/RSH-GT-5.7/0.63H 500UL SYRINGE	365J2441
007891	500F-RTC/RSH-GT-8.5/0.63H 500UL SYRINGE	365J2431
008150	1MF-RTC/RSH-GT-5.7/0.72LC 1ML SYRINGE	365K2811
008151	1MR-RTC/RSH-GT-5.7/0.72LC 1ML SYRINGE	365K2811
008155	1MF-RTC/RSH-EGT-HS-6.5/0.63H 1ML HS SYRINGE	365K2871
008655	2.5MF-RTC/RSH-EGT-HS-6.5/0.63H 2.5ML HS SYRINGE	365L2321
008950	10MF-RTC/RSH-GT-5.7/1.07LC 10ML SYRINGE	365N2721
008951	10MR-RTC/RSH-GT-5.7/1.07LC 10ML SYRINGE	365N2721

Information and support

Visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.