

Injection port | **General purpose blue**

Give your GC the seal of approval

Desired septa attributes:

- Reliably seal against the carrier gas pressure in the inlet.
- Capable of being pierced and resealing time after time without pieces of the septa being deposited in the GC inlet system.
- Does not contaminate or bleed material into the chromatographic system.



Features and benefits

- Combines significantly longer injection life, low bleed and low injection port adhesion.
- Packaged in glass jars for purity.

Recommended applications

Low temperature.

Product specifications

Sizes available (mm): 5, 6.4, 8, 9, 9.5, 10, 11, 11.5, 12.7, 17, and Shimadzu plug style.

Material	Durability	Resealing	Solvent resistance	Tear resistance	Maximum injector temperature
Silicone	Good	Good	Excellent	Good	275°C

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com.