

## Injection port | **AG3 low-bleed**

Give your GC the seal of approval

Desired septa attributes:

- Reliably seal against the carrier gas pressure in the inlet.
- Capable of being pierced and resealing time after time without pieces of the septa being deposited in the GC inlet system.
- Does not contaminate or bleed material into the chromatographic system.



### Features and benefits

- Combines significantly longer injection life, low bleed and low injection port adhesion.
- Packaged in blister packs (PK24,25) and glass jars (PK50,100) for purity.

### Recommended applications

Premium septa for autosamplers in all critical analyses.

### Product specifications

Sizes available (mm): 5, 6, 6.4, 7, 8, 9, 9.5, 10, 11, 11.5, 12.7, 17, and Shimadzu Plug style.

Material	Durability	Resealing	Solvent resistance	Tear resistance	Maximum injector temperature
Low-bleed silicone	Excellent	Excellent	Excellent	Excellent	300°C for 17 mm size 350°C for all other sizes

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com).