

Syringes | Headspace

Higher limits in headspace analysis

Part number	Volume	Needle length (mm)	Needle gauge	Needle OD (mm)	Needle ID (mm)	Needle tip	Pack size	Replacement plunger part number
CTC Analytics PAL3 RTC, Shimadzu AOC-6000, Thermo Scientific TriPlus RSH and Agilent PAL3								
Fixed needle syringes								
008155	1 mL	65	23	0.63	0.15	Side hole	1	032841
008655	2.5 mL	65	23	0.63	0.15	Side hole	1	032846
008855	5 mL	65	23	0.63	0.15	Side hole	1	032851
CTC Analytics and Thermo Scientific TriPlus								
Fixed needle syringes								
008131	1 mL	56	23	0.63	0.15	Side hole	1	031841
008136	1 mL	56	26	0.47	0.15	Side hole	1	031841
008631	2.5 mL	56	23	0.63	0.15	Side hole	1	031851
008636	2.5 mL	56	26	0.47	0.15	Side hole	1	031851
008831	5 mL	56	23	0.63	0.15	Side hole	1	031861

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.