

Injection port | **ConnectTite inlet liner**

Choose confidence in your analysis

Part number	Part description and detail
Agilent Intuvo 9000, 8890, 8860, 7890, 7890B, 7820A, 6890, 6850, 5890 and 4890; Bruker/Varian 1177; Lucidity GC-FID; PerkinElmer GC 2400, Clarus 690 and 590; Scion Instruments 8500 and 8300; Thermo Scientific TRACE 1600, 1610, 1300 series and GC Ultra	
092324	SGE ConnectTite inlet liners, length 78.5 mm PK5
092325	SGE ConnectTite inlet liners with top hole, length 78.5 mm PK5
092326	SGE ConnectTite inlet liners with bottom hole, length 78.5 mm PK5
0926232401	SGE OptChem™ ConnectTite inlet liner with pre-fitted CRS ONE o-ring, length 78.5 mm PK1
09262324	SGE OptChem ConnectTite inlet liners with pre-fitted CRS ONE o-ring, length 78.5 mm PK5
0926232501	SGE OptChem ConnectTite inlet liner with top hole and pre-fitted CRS ONE o-ring, length 78.5 mm PK1
09262325	SGE OptChem ConnectTite inlet liners with top hole and pre-fitted CRS ONE o-ring, length 78.5 mm PK5
0926232601	SGE OptChem ConnectTite inlet liner with bottom hole and pre-fitted CRS ONE o-ring, length 78.5 mm PK1
09262327	SGE OptChem ConnectTite inlet liners with bottom hole and pre-fitted CRS ONE o-ring, length 78.5 mm PK5
Bruker/Varian 1078, 1079, 1093 and 1094	
092030	SGE ConnectTite Bruker/Varian inlet liners 0.8 mm ID, length 54 mm (rest 0.25 mm) PK5
09203401	SGE ConnectTite Bruker/Varian inlet liner 0.8 mm ID, length 54 mm (rest 0.5 mm) PK1
092034	SGE ConnectTite Bruker/Varian inlet liners 0.8 mm ID, length 54 mm (rest 0.5 mm) PK5
Shimadzu Nexis GC-2030, GC-2025, GC-2014, GC-2010 Pro, GC-2010 and GC-17A	
092329	SGE ConnectTite Shimadzu inlet liners, length 95 mm PK5
092330	SGE ConnectTite Shimadzu inlet liners with top hole, length 95 mm PK5
092331	SGE ConnectTite Shimadzu inlet liners with bottom hole, length 95 mm PK5
Thermo Scientific TRACE GC Ultra, FOCUS GC and 8000 series	
092334	SGE ConnectTite Thermo Scientific inlet liners, length 105 mm PK5
092336	SGE ConnectTite Thermo Scientific inlet liners with bottom hole, length 105 mm PK5

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Agilent is a registered trademark of Agilent Technologies, Inc. Bruker is a registered trademark of Bruker Corporation. Shimadzu is a registered trademark of Shimadzu Corporation. Thermo Scientific is a registered trademark of Thermo Fisher Scientific, Inc.