

Injection Port | Inlet Liner Double Taper

Spot the Difference in the SGE Lineup!

Part number	Part description and detail
Agilent and Thermo Scientific TRACE 1300 Series GC	
09201801	Inlet Liner 4 mm ID Double Taper, Length 78.5 mm PK1
092018	Inlet Liners 4 mm ID Double Taper, Length 78.5 mm PK5
092230	Inlet Liners 4 mm ID Double Taper, Length 78.5 mm PK25

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com