

Injection port | Double taper inlet liner

Choose confidence in your analysis

Part number	Part description and detail
	Agilent Intuvo 9000, 8890, 8860, 7890, 7890B, 7820A, 6890, 6850, 5890 and 4890; Bruker/Varian 1177; Lucidity GC-FID; PerkinElmer GC 2400, Clarus 690 and 590; Scion Instruments 8500 and 8300; Thermo Scientific TRACE 1600, 1610, 1300 series and GC Ultra
09201801	SGE double taper inlet liner 4 mm ID, length 78.5 mm PK1
092018	SGE double taper inlet liners 4 mm ID, length 78.5 mm PK5
092230	SGE double taper inlet liners 4 mm ID, length 78.5 mm PK25
0926201801	SGE OptChem™ double taper inlet liner with pre-fitted CRS ONE o-ring, 4 mm ID, length 78.5 mm PK1
09262018	SGE OptChem double taper inlet liners with pre-fitted CRS ONE o-ring, 4 mm ID, length 78.5 mm PK5

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered trademark of Thermo Fisher Scientific, Inc.