

Separations | HPLC column ProteCol® C8

Turn imagination into results

Part number	Part description and detail
C8H	
120 Å Pore Size, Stainless Steel	
2C83-03L02SG	C8H x 3 µm Particle Size x 120 Å Pore Size, 50 mm Length x 2.1 mm ID, Stainless Steel with Guard ProteCol HPLC Column
2C83-03M46SG	C8H x 3 µm Particle Size x 120 Å Pore Size, 100 mm Length x 4.6 mm ID, Stainless Steel with Guard ProteCol HPLC Column
2C83-05N46S	C8H x 5 µm Particle Size x 120 Å Pore Size, 150 mm Length x 4.6 mm ID, Stainless Steel ProteCol HPLC Column
2C83-05N46SG	C8H x 5 µm Particle Size x 120 Å Pore Size, 150 mm Length x 4.6 mm ID, Stainless Steel with Guard ProteCol HPLC Column
2C83-05G02S	C8Hx 5 µm Particle Size x 120 Å Pore Size, 10 mm Length x 2.1 mm ID, Stainless Steel PK3 ProteCol HPLC Guard Column PK3
120 Å Pore Size - PEEK Coated	
2C83-03N02P	C8H x 3 µm Particle Size x 120 Å Pore Size, 150 mm Length x 2.1 mm ID, PEEK-Coated ProteCol HPLC Column
2C83-05N46P	C8H x 5 µm Particle Size x 120 Å Pore Size, 150 mm Length x 4.6 mm ID, PEEK-Coated ProteCol HPLC Column
2C83-05R46P	C8H x 5 µm Particle Size x 120 Å Pore Size, 250 mm Length x 4.6 mm ID, PEEK-Coated ProteCol HPLC Column
1000 Å Pore Size - Stainless Steel	
2C82-03N46S	C8H x 3 µm Particle Size x 1000 Å Pore Size, 150 mm Length x 4.6 mm ID, Stainless Steel ProteCol HPLC Column
1000 Å Pore Size - PEEK Coated	
2C82-03M02P	C8H x 3 µm Particle Size x 1000 Å Pore Size, 100 mm Length x 2.1 mm ID, PEEK-Coated ProteCol HPLC Column
2C82-03N02P	C8H x 3 µm Particle Size x 1000 Å Pore Size, 150 mm Length x 2.1 mm ID, PEEK-Coated ProteCol HPLC Column
2C82-03N46P	C8H x 3 µm Particle Size x 1000 Å Pore Size, 150 mm Length x 4.6 mm ID, PEEK-Coated ProteCol HPLC Column
1000 Å Pore Size - PEEKsil	
2C82-03T15K	C8H x 3 µm Particle Size x 1000 Å Pore Size, 10 mm Length x 150 µm ID, PEEKsil ProteCol HPLC Trap Column
2C82-03T30K	C8H x 3 µm Particle Size x 1000 Å Pore Size, 10 mm Length x 300 µm ID, PEEKsil ProteCol HPLC Trap Column

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com