

Separations | HPLC column ProteCol[®] Amino

Turn imagination into results

Part number	Part description and detail
120 Å pore size, stainless steel	
2AM-05R46S	Amino x 5 µm particle size x 120 Å pore size, 250 mm length x 4.6 mm ID, stainless steel ProteCol HPLC column
2AM-05R46SG	Amino x 5 µm particle size x 120 Å pore size, 250 mm length x 4.6 mm ID, stainless steel with Guard ProteCol HPLC column
120 Å pore size, PEEK coated	
2AM-05R46PG	Amino x 5 µm particle size x 120 Å pore size, 250 mm length x 4.6 mm ID, PEEK coated with Guard ProteCol HPLC column

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Separations | HPLC column ProteCol® Amino continued