

## Separations | GC column BP21 (FFAP)

Part number	Part description and detail
054462	0.22 mm ID x 0.25 µm film thickness x 25 m length BP21 GC capillary column
054465	0.25 mm ID x 0.25 µm film thickness x 30 m length BP21 GC capillary column
054466	0.25 mm ID x 0.25 µm film thickness x 60 m length BP21 GC capillary column
054468	0.32 mm ID x 0.25 µm film thickness x 25 m length BP21 GC capillary column
054471	0.32 mm ID x 0.25 µm film thickness x 30 m length BP21 GC capillary column
054469	0.32 mm ID x 0.25 µm film thickness x 50 m length BP21 GC capillary column
054472	0.32 mm ID x 0.25 µm film thickness x 60 m length BP21 GC capillary column
054480	0.32 mm ID x 0.5 µm film thickness x 50 m length BP21 GC capillary column
054474	0.53 mm ID x 0.5 µm film thickness x 25 m length BP21 GC capillary column
054477	0.53 mm ID x 0.5 µm film thickness x 30 m length BP21 GC capillary column
054475	0.53 mm ID x 1.0 µm film thickness x 15 m length BP21 GC capillary column
054478	0.53 mm ID x 1.0 µm film thickness x 30 m length BP21 GC capillary column

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)