

Model 513 Manipulation Couch

✓ V

✓ W

✓ Made

Overview:

A very strong and rigid couch designed for various manipulation and mobilisation procedures, this model features both positive and negatie adjustment of the foot and head sections with breathing hole and plug included as standard.



Benefits

- Lifetime frame & 5 year parts guarantee see website for further details
- · Proven lifespan in excess of 16 years
- Suitable for heavy manipulations & high usage
- Low running & maintenance costs
- Customisable to suit user's needs
- Safe working load of 250kgs (550lb / 39 stone)
 Optional foot bar switch system provides
 all-round, hands free height operation
- Extendable foot section option increase overall couch length for taller patients
- Head section can be lowered to -25° for head and neck procedures or for patient recovery

Overall Dimensions

Length: 200cm (79") approx

with foot extension 200 - 225cm (79" - 89") approx


Width: 63cm or 70cm

(25" or 28") approx

Height: 47cm - 97cm approx

(18" - 38") approx

Upholstery Dimensions



Options & Accessories:



Paper Roll Holder



Linked Castor Base



Foot Switch Operation



Foot Switch Extension



Dual Foot Switch Operation



Foot Bar Switch

NHS Approved Supplier





t: +44 (0)1449 767887 e: sales@plinthmedical.com www.plinthmedical.com

Operational Parameters



Model 513E - Electric Version

Model 513H - Hydraulic Version

Standard Features

- · Designed & manufactured in the UK
- Lifetime frame & 5 year parts guarantee
- · Very strong & rigid frame
- High-performance bearings on all pivot points
- Safe working load & lift capacity 250kg
- Hydraulic or electric height operation
- Retractable wheel system
- · Adjustable foot for uneven floors
- Anti-entrapment gas struts
- Head section adjustable from -25° to 80°
- Foot section adjustable from -20° to 85°
- Option of foot extension for taller patients
- Breathing hole and plug
- Choice of width 63cm or 70cm
- Hoist clearance 13cm

Technical Standards

- IP Rating IPx4
- EEC Classification MDD93/42/EEC Annex IX Class 1
- Electrical Standard IEC601 Part 1 + 60601 P (EN)
- Vinyl Flame Retardant to CRIB 5 (BS 7176)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]
- CE Marked
- Externally audited to ISO9001:2015 by LRQA

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service may be available if required
- Overseas delivery 7-10 days plus shipping time

Upholstery Vinyl

- Flame retardant to BS 5852 Ignition Source 5 CRIB 5 (BS 7176) Medium Hazard
- Independently tested to verify effectiveness at reducing the spread of MRSA, E.coli & other similar bacteria & fungi
- Rated at 500,000 Martindale rubs without damaging or marking the top layer, making it suitable for all high-traffic, heavy-duty areas
- · Phthalate free
- Please note: If the upholstery colour is not specified, 'Almond' will be supplied



Plinth Medical Ltd Barric Lane Occold, Eye Suffolk, United Kingdom IP23 7PX t: +44 (0)1449767887

e: sales@plinthmedical.com

www.plinthmedical.com

