

# RETENTION+™

## SCULPTING LIQUID

### WHAT IT IS

A keratin-bonding, cross-linking sculpting liquid for lift-prone clients.

### WHAT IT DOES

Creates enhancements with superior adhesion to the natural nail.

### WHY YOU NEED IT

- Use on clients with enhancements that are prone to lifting.
- Excellent adhesion makes primer use optional.

### DIRECTIONS FOR USE

- Perform P.R.E.P. and ensure the nail is thoroughly cleansed and prepared.
- Apply NailPrime, if desired.
- Sculpt the nail in 3 zones using a mix ratio of 1.5 parts liquid to 1 part powder.
- Shape and refine using **CND™** files and buffers.
- Visit [cnd.com](http://cnd.com) for complete Step-by-Step instructions.



### AVAILABLE SIZES

- Professional: 118 mL (4 fl oz)  
236 mL (8 fl oz)  
473 mL (16 fl oz)  
946 mL (32 fl oz)  
1894 mL (64 fl oz)  
3785 mL (Gallon)