



PARADISE
DENTAL TECHNOLOGIES

(406) 626-4153

customer.care@pdt dental.com

8275 US Highway 10 W Missoula, MT 59808

PINEYRO ARCH™ KIT

Created in collaboration
with Dr. Alfonso Pineyro



EXCELLENT HYGIENE ACCESS FOR FULL FIXED IMPLANT RESTORATIONS

Your R930 KIT includes
a T005A-F **FlipTop™ Cassette**
and four instruments:

- R630 **PINEYRO ARCH TI 1**
Anterior Arch
- R631 **PINEYRO ARCH TI 2**
Posterior Arch
- R632 **PINEYRO ARCH TI 3**
Mesial and Distal
- R633 **PINEYRO ARCH TI 4**
Ideal abutment cleaning

*No need for
resharpening.*

Design follows
the circular
anatomy of dental
implants

Greater
efficiency
and less
gingival trauma

No need for
removal of fixed
restorations



Excellent hygiene access for
full fixed implant restorations.
with no need for resharping.

-Dr. Alfonso Pineyro






CELLERANT
BEST OF CLASS
HYGIENE AWARD

2021

Dental **PRODUCT
SHOPPER**

**BEST
PRODUCT**



Pineyro Arch Ti 1

Anterior Arch - to debride the anterior placed implants for a full-fixed hybrid, especially lingual area of lower anterior teeth, where calculus builds up most frequently




Pineyro Arch Ti 2

Posterior Arch - to debride the posterior placed implants for a full-fixed hybrid, the design of the shank allows for proper access



Pineyro Arch Ti 3

Mesial/Distal Arch - to debride the posterior hard-to-reach mesial/distal areas on the implants for a full-fixed hybrid



Pineyro Arch Ti 4

Specialty Arch - or all implants prior to taking impressions - clean off abutments for single to hybrid cases prior to cementation



Go here for more information or to watch a demo video by Dr. Pineyro.