TECHNICAL DATA SHEET



SybraChem Alcohol Sanitiser Liquid 70%

A quick-drying sanitiser gel that provides immediate sterilisation without leaving a sticky residue on your hands. Use this product to supplement routine hand washing and drying. This product is ideal for use before handling food or working with patients, where absolute cleanliness is essential.

PROPERTIES

- Has no residual fragrance or harmful residue.
- Leaves hands feeling smooth and silky.
- Not sticky and non-greasy.

APPLICATION

Rub the hand sanitiser gel on hands while spreading evenly on the skin in a washing motion for approximately 1 minute to allow the formulation time to kill surface bacteria and evaporate naturally. Do not wipe off.

TECHNICAL INFORMATION

Product Name: SybraChem Alcohol Sanitiser Liquid 70%		
TEST	UNIT	SPECIFICATION
Appearance		Clear liquid
Density	g/cm ³	0.90 – 0.94
рН		8.00 – 9.00

HEALTH, SAFETY, AND THE ENVIRONMENT: For full **Health and Safety** guidance, obtain a **Material Safety Data Sheet** from your representative, or online at www.sybrachem.africa. The recommendations contained herein are based on infield use experience and are, to the best of our knowledge, accurate. Since terms and conditions of actual use are beyond our control, all recommendations are made without warranty, express or implied.

DISCLAIMER: SybraChem (Pty) Ltd believes that this information is correct at the date of publication but does not warrant the accuracy thereof. The use of the product designated herein will be at the sole risk of the user. SybraChem (Pty) Ltd will not accept any liability for any loss or damage, howsoever caused by or arising from the use of this information or the use, application, adaptation, or processing of the product described herein or its use in combination with any product, material or in any process.