

LEVEL 2 SURGICAL FACE MASK

- TGA LISTED (50 PACK)



This Level 2 TGA Listed Surgical Face Mask has superior Bacterial Filtration Efficiency (BFE) of >98% and is fluid resistant up to 120mmHg. This mask is a suitable face mask for medical professionals who may have exposure to body and body fluids as well as risk of potential infection. This mask can also be used in industries such as cleaning and maintenance. With 3 layers of protection to filter harmful substances and provide protection against biological aerosols, dust, germs, haze, allergens, smoke, spores, microorganisms, pollution, ash, pollen, garden materials, and odours.

- Colour: White
- Size: one size fits all
- Complies with ASTM Level 2
- CE & FDA Approved
- Material: 3-ply non-woven fabric
- Product Dimensions: 17.5 x 9.5 cm
- Fluid resistant
- Hypoallergenic
- >98% BFE
- TGA Listed
- Single use only

Features ASTM Level 2:

High Fluid Resistance	120 mmHg
Filtration Efficiency	BFE ≥ 98%
Breathability - Delta P	0.1 micron < 5.0 mm H ₂ O/cm ²
Flame Spread	Class 1

Time Use Limitation: Single use only.

Storage: Store in a well-ventilated place. Protect from Sunlight.

FITTING INSTRUCTIONS



1

Hold mask in your hands as shown, with the nosepiece at your fingertips and the headbands hanging freely below your hand.




2

Place mask under your chin with the nosepiece up. Pull the ear loops behind your ears.

Caution:

The mask should be fit checked for each wearing. Place both hands completely over the mask and exhale. If air leaks around nose or the edges, adjust the nosepiece and/or earloops until a good fit is achieved.



**Outer layer -
Non-Woven material**
(repels fluid and bacteria)

**Middle layer -
Meltblown
Non-Woven**
mould filter to absorb
fine particles in the air

**Inner layer -
Non-Woven**
protection layer for
double filtration