

NeoLok Faraday Bag with Charging Kit

The Mission Darkness™ NeoLok Faraday Bag with Charging Kit is designed to keep devices live during evidence transport and up to the point of data extraction. Maintaining power to the device may allow for recovery of more data, and faster brute force unlock times. The shielded compartment will hold a device along with the included power bank and cables, while the pocket on the back stores accessories. The transparent window located on the front of the bag allows the user to see the power bank LED indicators showing power bank charge status.

Kit Includes

- NeoLok Magnetic Closure Faraday Bag
 - Portable Power Bank (capacity may vary depending on kit configuration)
 - USB Extension Cable
 - Lighting Cable
 - USB-C Cable
-

Usage Instructions

1. Connect a device to the USB power bank using the correct charging cable
 2. Place both the device and power bank into the main faraday bag shielded compartment, with the battery indicator light facing forward (towards the window). The battery indicator should stay lit, but if it does not, press the top of the battery through the bag to activate the indicator. It's best to check the location of this button before using the kit
 3. Close the bag by folding the top of bag over twice to pair the magnets together
 4. Use the button snap to ensure the bag will stay closed and shielded during transport
 5. Place the evidence card in the back clear pocket
-



Tel: (805) 318-3212
info@mosequipment.com
www.mosequipment.com

MADE IN USA