

- EZBio® Carboys are Gamma Sterilized and double bagged in a Cleanroom environment
- VersaCap® technology provides the flexibility of interchangeable connectors
- EZBio® Carboy body features advanced ergonomic design
- EZfil® wide-mouth design for easy filling and cleaning

Caution

- Do not use carboys under pressure or full vacuum
- Do not use carboys at temperatures greater than the rated maximum
- Do not store oxidizing agents in carboys
- Do not flame carboys
- Do not expose to organic compounds. Carboy will become prone to stress fracturing

Made in the USA



ISO 13485


Foxx Life Sciences' EZBio® Carboy is an advanced fluid handling system developed with the end-user in mind. The EZBio® Carboy serves as a highly versatile container, suitable for all fluid management applications. Carboys are made entirely of USP Class VI, FDA 21 CFR 177 Grade materials and are available Sterile or Non-Sterile, double bagged. All EZBio® Carboys are carefully triple inspected and double bagged in an ISO 7 cleanroom, in a facility certified to ISO 13485:2016 device manufacturing.

These carboy come in a distinctive rectangular shape that saves valuable bench space and features a wide-mouth design for easy filling and cleaning. The carboy body is ergonomically designed, with large top handles and molded hand grips in the base, allowing the user to easily hold and manipulate during filling, lifting, and pouring. Additionally, the EZBio® Carboy can be coupled with our patented VersaCap® Technology, providing interchangeable and leak-proof top connections. Our caps are molded not machined, to eliminate particles or weakness. Choose from a wide variety of adapter inserts to match your application needs. The EZBio® Carboy features easy-to-read metric graduation marks certified to ±5% accuracy. The carboy material is molded into the side for easy identification, and tells the user if it can be Autoclaved or not.

The Polycarbonate Carboy is made from a clear USP Class VI, FDA Grade autoclavable material. This Carboy is suitable for aqueous solutions, as well as most life sciences buffer solutions. The recommended autoclaving cycle is 15 PSI for 20 minutes at 121°C.



- 1 Novel shape maximizes bench and shelf space
- 2 Bottom handles provide a secure grip while lifting and transporting
- 3 Molded-in carboy identification marks
- 4 VersaCap® technology with versatile adapters and replaceable gaskets
- 5 EZBio® body design allows for maximum solution recovery
- 6 EZfil® wide-mouth shape simplifies filling and cleaning



Product Number Sterile	Product Number Non-Sterile	Description	Capacity (L)	Material	Cap Size (mm)	Qty/Pack
518-5C51-FLS	518-5C52-FLS	Carboy, PC, 2.5L, VersaCap® 83B	2.5L	PC	83B	1
518-5F51-FLS	518-5F52-FLS	Carboy, PC, 5L, VersaCap® 83B	5L	PC	83B	1
518-5G51-FLS	518-5G52-FLS	Carboy, PC, 10L, VersaCap® 83B	10L	PC	83B	1
518-5H61-FLS	518-5H62-FLS	Carboy, PC, 20L, VersaCap® 120mm	20L	PC	120mm	1
518-5J61-FLS	518-5J62-FLS	Carboy, PC, 40L, VersaCap® 120mm	40L	PC	120mm	1
518-9G52-FLS	518-9G51-FLS	Round Carboy, PC, 10L, VersaCap® 83B	10L	PC	83B	1
518-9H52-FLS	518-9H51-FLS	Round Carboy, PC, 20L, VersaCap® 83B	20L	PC	83B	1



BPA-Free



EZBio® Carboys

Foxx carries the largest range of carboys in the world, offering versatility in both material and volume. EZBio® Carboys come in your choice of High-Density Polyethylene (HDPE), Polypropylene (PP), Polycarbonate (PC), and Polyethylene Terephthalate Glycol (PETG). Sizes range from 2.5L (0.5 Gal) to 75L (20 Gal). Amber and Dark Amber options are also available for HDPE and PP.