




- Sold in 15ml and 50ml centrifuge tube sizes
- 0.1µm, 0.2µm, and 0.45µm PES membrane pore-size options
- Permanent vacuum hose connection to minimize set-up time
- Made from Class VI Materials
- Manufactured in a Class 100,000 Clean Room
- Manufactured according to ISO 13485 standards in North America
- **Gamma Sterilized in the USA** 

Foxx Life Sciences has developed the world's first hands-free centrifuge tube filtration system. Each filtration unit includes a sterile centrifuge tube and cap, 50ml funnel with PES membrane, funnel adapter (for the 15ml tube assembly), vacuum hose adapter, and dual centrifuge tube stand.

The **Autofil®** Centrifuge Funnel system is used for sterilization of tissue cultures and microbiological media, buffers, reagents, and other laboratory solutions. This system avoids potential contamination by eliminating any unnecessary fluid transfer to the centrifuge tube post-sterilization, and reduces risk of spills.

This hands-free filtration system works in conjunction with the Quick-Connect Pedestal or Cradle Ring which provides on/off vacuum flow control. These two filtering options provide a permanent vacuum hose connection if a Quick-Connect Pedestal or Cradle Ring is not available.


The **Autofil®** Centrifuge Funnel system is offered as a full assembly, or funnel only. Available sizes include a 15ml and 50ml centrifuge tube with 0.1 µm, 0.2 µm, and 0.45 µm PES membrane pore-size options.



15ml **Autofil®** Centrifuge Funnel system with Quick-Connect Pedestal base (sold separately)

ISO 9001
ISO 13485

- 1 50ml filtration funnel with forward-facing and raised graduation marks for easy reading
- 2 Available in 0.1µm, 0.2µm, and 0.45µm asymmetric membrane filter pore sizes
- 3 15ml **Autofil®** Centrifuge Funnel system
- 4 50ml **Autofil®** Centrifuge Funnel system
- 5 Vacuum hose adapter
- 6 Quick-Connect pedestal base for hands-free filtration





	Item Number	Volume (ml)	Membrane	Filter Diameter (mm)	Description
Full Assembly	146-5113-RLS	50	0.1 µm PES	40	50 mL Centrifuge Tube and Funnel with 0.1 µm PES Membrane, 24/pk
	1402-RLS	50	0.2 µm PES	40	50 mL Centrifuge Tube and Funnel with 0.2 µm PES Membrane, 24/pk
	1403-RLS	50	0.45 µm PES	40	50 mL Centrifuge Tube and Funnel with 0.45 µm PES Membrane, 24/pk
	146-4113-RLS	15	0.1 µm PES	40	15 mL Centrifuge Tube and Funnel with 0.1 µm PES Membrane, 24/pk
	1411-RLS	15	0.2 µm PES	40	15 mL Centrifuge Tube and Funnel with 0.2 µm PES Membrane, 24/pk
	1412-RLS	15	0.45 µm PES	40	15 mL Centrifuge Tube and Funnel with 0.45 µm PES Membrane, 24/pk
Funnel-Only	146-2113-RLS	50	0.1 µm PES	40	50 mL Funnel Only with 0.1 µm PES Membrane, 48/pk
	146-2213-RLS	50	0.2 µm PES	40	50 mL Funnel Only with 0.2 µm PES Membrane, 48/pk
	146-2313-RLS	50	0.45 µm PES	40	50 mL Funnel Only with 0.45 µm PES Membrane, 48/pk
	146-1113-RLS	15	0.1 µm PES	40	15 mL Funnel Only with 0.1 µm PES Membrane, 48/pk
	146-1213-RLS	15	0.2 µm PES	40	15 mL Funnel Only with 0.2 µm PES Membrane, 48/pk
	146-1313-RLS	15	0.45 µm PES	40	15 mL Funnel Only with 0.45 µm PES Membrane, 48/pk
Hard-ware	1116-RLS				Vacuum Filtration Pedestal
	1122-RLS				Vacuum Filtration Cradle Ring
Acces-sories	146-9004-RLS				Polystyrene Dual Centrifuge Tube Stand, 5/pk

The Most Patented Vacuum Filtration System in the World!

US Patent Numbers: 7,798,333, 8,158,009, 8,235,221, 8,231,012, 8,157,104, 8,480,895, 7,011,755, 7,661,538, D587580
 Also Patented in Europe, India, China, and Japan