

**RAW ORGANIC WHEY
RAW ORGANIC WHEY
BOX 254
FAIRFIELD IA 52556**


REPORT OF ANALYSIS

For: (34261) RAW ORGANIC WHEY
Raw Organic Whey
Raw Grass Fed Whey
RAW01A

Analysis	Level Found	Units	Reporting		Analyst- Date	Verified- Date
	As Received		Limit	Method		
Sample ID: 1A	Lab Number: 13285383					
Arsenic (total)	n.d.	ppm	0.10	USP <233>(ICP-MS)	trh1-2019/08/27	kkh9-2019/08/27
Cadmium (total)	n.d.	ppm	0.020	USP <233>(ICP-MS)	trh1-2019/08/27	kkh9-2019/08/27
Lead (total)	n.d.	ppm	0.10	USP <233>(ICP-MS)	trh1-2019/08/27	kkh9-2019/08/27
Mercury (total)	n.d.	ppm	0.01	USP <233>(ICP-MS)	trh1-2019/08/27	kkh9-2019/08/27

All results are reported on an AS RECEIVED basis., n.d. = not detected , ppm = parts per million, ppm = mg/kg

For questions please contact:



Kailey Parr
Account Manager
kparr@midwestlabs.com (402)829-9863

The result(s) issued on this report only reflect the analysis of the sample(s) submitted.

Our reports and letters are for the exclusive and confidential use of our clients and may not be reproduced in whole or in part, nor may any reference be made to the work, the results, or the company in any advertising, news release, or other public announcements without obtaining our prior written authorization.

REPORT NUMBER

19-239-9488

REPORT DATE
Aug 27, 2019
RECEIVED DATE
Aug 21, 2019

SEND TO
34261



PAGE 2/2
ISSUE DATE
Aug 27, 2019

**RAW ORGANIC WHEY
RAW ORGANIC WHEY
BOX 254
FAIRFIELD IA 52556**

REPORT OF ANALYSIS

For: (34261) RAW ORGANIC WHEY
Raw Organic Whey
Raw Grass Fed Whey
RAW01A

Detailed Method Description(s)

ME 081

Sample analysis is conducted by ICP-MS which follows an acid digestion/preparation of the sample which destroys and solublizes the sample. The ICP-MS analysis uses a plasma to induce energy into prepared samples so as to breakdown the compounds present and create a stream of elemental ions. The ions are then separated by a mass spectrometer in to their individual elements. The mass spectrometer measures the masses of the elements present and quantifies the levels present. These results are correlated to known levels of standards and calculated back to original concentration in the sample analyzed.

ME 080

Sample preparation for metals analysis referenced by <USP 233> follows MWL ME 080 which is a microwave assisted wet-ash digest or a "neat" dilution.

The result(s) issued on this report only reflect the analysis of the sample(s) submitted.

Our reports and letters are for the exclusive and confidential use of our clients and may not be reproduced in whole or in part, nor may any reference be made to the work, the results, or the company in any advertising, news release, or other public announcements without obtaining our prior written authorization.