



4300 Series

LOOPED END MOPS

PRODUCT DATA SHEET

TYPICAL APPLICATIONS: HEALTHCARE



DESCRIPTION

The Sterling Plus series is constructed of silver treated cotton/rayon/synthetic fibers. Sterling Plus antimicrobial loop mops would be ideal for the commercial laundry rental market, for food service, healthcare, and for all Jan/San applications. Cotton/rayon/synthetic fibers are treated with a silver antimicrobial agent which will not wash out into the environment. The antimicrobial agent protects the mop head from bacterial and fungal growth, which inhibits yarn degradation and offensive odors.

PRODUCT SPECS

O'Dell Item	Size	Weight (lbs)	Headband	Tailband Spread	Working Mop Length	Total Mop Length
4300S	Small	0.75	Yellow 5" Mesh	15"	15"	30"
4300M	Medium	1.00	Green 5" Mesh	16"	17"	34"
4300L	Large	1.35	Red 5" Mesh	18"	21"	42"
4300XL	Extra Large	1.70	Blue 5" Mesh	22"	23"	46"
Color Options:	Blue, Green, Orange, White		Launderability:	Launderable		

PACKAGING SPECS

O'Dell Item	Case Pack Qty	Case Weight (lbs)	Case Dimensions	Case Cube
4300S	12	10.0	19.75 x 9.25 x 9.25	0.98
4300M	12	13.5	19.75 x 11.63 x 11.63	1.55
4300L	12	18.0	19.75 x 11.63 x 11.63	1.55
4300XL	12	21.0	19.75 x 11.63 x 14.5	1.93

Note: Each mop is individually polybagged

DIRECTIONS FOR USE

Remove mop from polybag and attach to either a jaws or quick change mop handle. Submerge mop into prepared cleaning solution to fully saturate. Wring mop to desired level of saturation and begin mopping using a figure eight motion. If mop becomes too dry or soiled before completing the area, dispose of old solution and replace with clean solution.

PRODUCT CARE

Thoroughly rinse. Avoid long term exposure to harsh chemicals and abrasive surfaces. After every use wash and ring out mop. Store mops by hanging in a warm, dry ventilated room.

ADDITIONAL OPTIONS

Narrow Band, Industrial Flat, Twist On, Scrub Pad, Cut End