




# 5500 Series

## LOOPED END MOPS

## PRODUCT DATA SHEET

**TYPICAL APPLICATIONS: GENERAL MAINTENANCE AND HEALTHCARE**



### DESCRIPTION

4-ply, premium quality, synthetic blend mophead for healthcare applications. Prelaundered with limited break in.

### PRODUCT SPECS

O'Dell Item	Size	Weight (lbs)	Headband	Tailband Spread	Working Mop Length	Total Mop Length
5500S	Small	0.75	Yellow 5" Mesh	15"	12.5"	25"
5500M	Medium	1.00	Green 5" Mesh	16"	14.5"	29"
5500L	Large	1.35	Red 5" Mesh	18"	18.5"	37"
5500XL	Extra Large	1.70	Blue 5" Mesh	22"	20.5"	41"
<b>Color Options:</b>	Blue, Green, White, Red, Yellow, Orange		<b>Launderability:</b>	Launderable		

### PACKAGING SPECS

O'Dell Item	Case Pack Qty	Case Weight (lbs)	Case Dimensions	Case Cube
5500S	12	10.0	19.75 x 9.25 x 9.25	0.98
5500M	12	13.5	19.75 x 11.63 x 11.63	1.55
5500L	12	18.0	19.75 x 11.63 x 11.63	1.55
5500XL	12	21.0	19.75 x 11.63 x 14.5	1.93

Note: Each mop is individually polybagged

### DIRECTIONS FOR USE

Remove mop from polybag and attach to either a jaws or quick change mop handle. Submerge mop into prepared cleaning solution to fully saturate. Wring mop to desired level of saturation and begin mopping using a figure eight motion. If mop becomes too dry or soiled before completing the area, dispose of old solution and replace with clean solution.

### PRODUCT CARE

Thoroughly rinse. Avoid long term exposure to harsh chemicals and abrasive surfaces. After every use wash and ring out mop. Store mops by hanging in a warm, dry ventilated room.

### ADDITIONAL OPTIONS

Narrow Band, Industrial Flat, Twist On, Scrub Pad, Cut End