



Or not to delete and block?

#73 – 5 Options to Address Toxicity in Social Media with Teryn Darling

Join me in the latest episode of Tatter-a-fact PMU Podcast where I dive into a pressing topic: toxicity on social media. I discuss the five approaches to handling this, from **constructive criticism to outright negativity**, weighing the **pros and cons** of each. Plus, learn about the **benefits of Liquidcaine numbing!** Tune in now.

These are topics I talk about:

- [Addressing the Toxic and Negative Aspect of Social Media](#)
- [Your 5 Options When Addressing Toxic Social Media Behavior](#)
- [How to Have a Dialogue with Controversial Comments](#)
- [All About Liquidcaine Numbing](#)

[Watch the Full Episode!](#)



Li Pigments Aqua - Hazelnut

Temp: Neutral

Ideal Fitzpatrick: I, II & III

Base: More neutral Yellow & Green base, contains a small amount of Red.

Type: Inorganic & Organic

An absolute MUST HAVE for all artists! Single needle artist? You'll LOVE this color for your blonde clients!! A fan of bigger needles? You'll still love this color and find it heals to an amazing.gorgeous light to medium brown depending on skin type and needle choice.

[Try Hazelnut Today!](#)

4 DIFFERENCES BETWEEN ORIGINAL AND GLOBAL AQUA® & VELVET® LINES

- 1 Slight changes in pigment appearance
- 2 Minor differences in pigment viscosity
- 3 Similar implantation & performance
- 4 Slight visual changes in healed results

PMU In the Know with Darlene and Mytia

PMU In the Know with Darlene and Mytia

Darlene, Li Pigments president and founder, established the company over 30 years ago. Her daughter, Mytia, serves as vice president and senior chemist. With their extensive expertise in pigments and PMU, they're excited to share valuable insights! Look out for condensed tips in every Gazette, with detailed blogs for further reading.

Are there any distinct differences between the Aqua® & Velvet® GLOBAL and the Aqua® & Velvet® ORIGINAL?

Pigment Appearance: While there is a slight change in the appearance of the GLOBAL shade in the bottle, Li Pigments have made every effort to maintain the base of the ORIGINAL finished pigment counterpart shade.

Pigment Viscosity: Some GLOBAL shades may exhibit minor differences in viscosity from the ORIGINAL counterpart shade. Depending on the pigment, one may experience a slightly thinner consistency.

Pigment Implantation & Procedure Performance: The new GLOBAL pigments have performed excellently in the field comparable to ORIGINAL. With the GLOBAL, the live model testing results determine impressive pigment implantation, crisp lines, strokes that maintain without blurring, strong retention, concentrated levels and lasting results.

Healed Results: With the GLOBAL, technicians will not observe remarkable visual changes in healed results. However, a slight shade shift in appearance may be noticed, with some colors appearing deeper or lighter in tone from the ORIGINAL counterpart shade.

[Read the full blog](#)

[Get Gray Vanish!](#)

The Ultimate Scalp Micropigmentation Course

2 DAY HANDS ON CLASS WITH MODELS

When:
Nov. 4th-5th, 2024
AND
Nov. 6th-7th, 2024

Where:
Girlz Ink Studio
Las Vegas, Nevada

SCALP SISTA'S
JESSICA ABBOTT & KERRY WELSTEAD

For further information email training@girlzink.com or call (702)809-8237

[SHOP THE STORE!](#)

[LEARN ON THE ACADEMY!](#)

