

Available in a variety of sizes, medical sheepskins are ideal for surgical and long-term patients who require buttock, shoulder, elbow, heel and/or lower back protection. Medical sheepskins unique qualities help prevent the three main causes of tissue breakdown: Pressure, Moisture, and Shearing Action.

The considerable cost burden to hospitals and nursing homes for the treatment of pressure ulcers is dramatically reduced with the use of medical sheepskins to prevent ulcers and assist in the recovery of patients with existing ulcers.

The real breakthrough has been the development of High Temperature U.R. (Urine Resistant) Medical Sheepskins. When cared for with SKINSAN, they are safely washable at temperatures of 80°C with no deterioration, thermally and chemically disinfecting them, allowing for their safe use many times.

These bed pads are 100% natural Australian Merino Sheepskin, available in medical green only. Wash at high temperature (80°C) using SKINSAN Liquid Sheepskin Detergent, and tumble dry.



Chair and Bed Pad
Approximately 39" by 30"

MSHE



Bed Pad

MSHP4 48" by 30"

MSHP5 60" by 30"



Bed Wrap
As a pillow cover, for shoulders, coccyx, or heels and ankles.

BWM1



Wheel Chair Pad
Contour Shape

MC1 21" by 18"



Wheel Chair Pad
Standard Shape

MC2 17" by 15"



Wheel Chair Covers
for seat and back

UNICW18M 16" by 18" (seat and back)
ARI Arm Rest (pair)-measurements required
CPI Calf Pad (pair)
FRI Foot Rest (pair)



Adjustable Lumbar
Support

LS3



Australian Sheepskin Apparel

provides sheepskins of the highest quality which meet Australian and international standards.

Products are available in both high quality medical and domestic grades.

Medical sheepskins may be laundered at high temperatures and tumble dried. Sheepskins may be washed repeatedly when laundered to the recommended standards, and can be sanitized to the same degree as normal hospital linen.

Please ask about Australian Sheepskin Apparel special sanitizing detergents.



SKINSAN
Detergent

#1. SSC0.2 2 oz. #3. SSC1 1 litre
#2. SSC2.5 250 ml #4. SSC4 4 litres
#5. SSC20L 20 litres
#6. HP2-1 Hand pump for bulk container

Warranty void if not laundered in SKINSAN Natural Sheepskin Detergent.


Australian Sheepskin Apparel
Head Office: Saskatoon Saskatchewan

Phone: (306) 934-7119 Fax: (306) 934-1597

Email: playzad@sasktel.net

www.australianmedicalsheepskins.com

Available from:



Printed in Canada



Australian Sheepskin Apparel

"Canada's Medical Sheepskin Experts"



Pressure relief, moisture absorption and shear reduction you can rely on.

The most naturally wonderful, comfortable and effective product for the prevention and treatment of wounds that health care has ever seen. Once you explore and feel the texture of Australian Medical Sheepskin, you will understand the amazing benefits to both patient and caregiver alike. This product was specifically designed to save health care millions of dollars by not only treating but preventing bed sores, while at the same time giving comfort to long-term and palliative care patients.



The billiard balls show the superior support capacity of Australian Medical Sheepskin as compared to an imitation.

Look for the Leather back...

Only real sheepskin has a leather back. Manufactured fabrics have a knit back that may or may not be coated with a rubber or plasticized material.

If you don't see it here we will endeavor to accommodate your custom order.

Soft velcro securing straps won't cut into skin and allow for easy size adjustment.



Limb Pressure Guard
(4 sizes)

LPG1 - 4" LPG3 - 6"
LPG2 - 5" LPG4 - 8"



Heel Protectors
(pair)

UNIHPM



Palm Protectors
(pair)

PPM - Regular
PPMF - with Finger Separators



Elbow Protectors
(1 each)

EPSM



Neck Support Ring
(2 sizes)

NRM3 (3 inches)
NRM4 (4 inches)




Crutch Covers
(pair)

UACP (under arm)
HGCP (hand grip)



Pressure Care Slipper
with Rubber Sole

PCBMHS (small) PCBMHL (large)
PCBMHM (med) PCBMHXL (X-large)



Surgical Slipper
(open toe)

UNISSMM
(one size fits all)



Reusable Coccyx Panties
and Pads

RCP-2 (small - medium)
RCP-3 (medium - large)




Coccyx Pads (5" by 5")

CP1-Pack of 6
CP2-Pack of 12



Wrap Around Slipper
with Rubber Sole

EPSMHS (small) EPSMHL (large)
EPSMHM (med) EPSMHXL (X-large)



Medical Moccasin
with suede or rubber sole

MMRS (rubber sole)
MMSS (suede sole)



Medical Mules
(4 sizes)

SSMS (small) SSML (large)
SSMM (med) SSMXL (X-large)



Executive Mules
for style and benefits of Medical Sheepskin

ESMS (small) ESML (large)
ESMM (med) ESMXL (X-large)



Short Boot
(easy access - 4 sizes)

SBMS (small) SBML (large)
SBMM (med) SBMXL (X-large)



Regular Long Boot
(easy access - 4 sizes)

LBMS (small) LBML (large)
LBMM (med) LBMXL (X-large)



Domestic Slippers

NC31 Available in mens and ladies sizes with suede or rubber soles.



Inner Soles

NC311S Available in mens and ladies sizes.




Short Split Boot
For Hammer Toe & Diabetic ulcers

SSBS (small) SSBL (large)
SSBM (med) SSBXL (X-large)



Long Split Boot
Adjusts to suit the leg & foot

LSBS (small) LSBL (large)
LSBM (med) LSBXL (X-large)



Baby Cozy Steps

BCSD



Domestic Sheepskin Rugs
Natural Shape (approx. 34" by 72")

White, Champagne, Beige,
BRD