

Available in a variety of sizes, medical sheepskins are ideal for surgical and long-term patients who require buttock, shoulder, elbow, heel and/or lower back protection. Medical sheepskins unique qualities help provide the four main requirements the human body needs to prevent wounds or treat pre existing wounds with moisture absorption, pressure reduction, shearing reduction and maintenance of standard body temperature 98.3 degrees.

The considerable cost burden to hospitals and nursing homes for the treatment of pressure ulcers is dramatically reduced with the use of medical sheepskins to prevent ulcers and assist in the recovery of patients with existing ulcers.

The real breakthrough has been the development of High Temperature U.R. (Urine Resistant) Medical Sheepskins. When cared for with SKINSAN, they are safely washable at temperatures of 80°C with no deterioration, thermally and chemically disinfecting them, allowing for their safe use many times.



Crutch Covers* (pair)

UACP (under arm)
HGCP (hand grip)



Walker Grips

WG1 – Walker Grips
* Supply Measurements



Knee Walker with elastic skirt 15" by 7"

KW1



Adjustable Lumbar Support

LS3



Domestic Long Hair Throw Rugs

DLH



Domestic Rug
Natural Shape (approximately 39" by 30")

Champagne colour.
BRD



Fire Retardant Sheepskin Hood

FRSH
SM, MED, LG, XL



Limb Pressure Guard (4 sizes)

LPG1 - 4" LPG3 - 6"
LPG2 - 5" LPG4 - 8"



Floor Rugs

FR4



Car Seat Cover

CSC



Boot Liners (pair)

WBL
SM, MED, LG, XL
mens/womens



Inner Soles

NC31IS Available in mens and ladies sizes.



Australian Sheepskin Apparel

provides sheepskins of the highest quality which meet Australian and international standards.

Products are available in both high quality medical and domestic grades.

Medical sheepskins may be laundered at high temperatures and tumble dried. Sheepskins may be washed repeatedly when laundered to the recommended standards, and can be sanitized to the same degree as normal hospital linen.

Please ask about Australian Sheepskin Apparel special sanitizing detergents.

Warranty void if not laundered in SKINSAN Natural Sheepskin Detergent.



SKINSAN Deergen

#1. SSC0.2 2 oz. #3. SSC1 1 litre
#2. SSC2.5 250 ml #4. SSC4 4 litres

*Crutch, Wheelchair Accessories and Specialized Footwear are GST exempt.


Australian Sheepskin Apparel

Head Office: Saskatoon Saskatchewan
Phone: (306) 934-7119 Fax: (306) 934-1597
Email: australiansheepskinapparel@shawbiz.ca

Please view our series of **SHORT** instructional videos detailing the proper use and benefits of using **Australian Medical Sheepskin** on our **YouTube Channel** through the link available on our website:

www.australianmedicalsheepskins.com

Available from:



Printed in Canada



"Canada's Medical Sheepskin Experts"



These bed pads are 100% natural Australian Merino Sheepskin, available in medical green only. Wash at high temperature (80°C) using SKINSAN Liquid Sheepskin Detergent, and tumble dry.



Chair and Bed Pad
Approximately 39" by 30"



MSHEP Bed Pad
Approx. 85" by 30"

MSHE



Bed Pad

MSHP4 48" by 30"

MSHP5 60" by 30"



Bariatric Bed Pad
Approximately 80" by 60"

BBP - Bariatric Overlay



Wheel Chair Covers*
for seat and back

UNICW18M - 16" by 18" (seat and back)
ARI - Arm Rest (pair)-measurements required
CPI - Calf Pad (pair) FRI - Foot Rest (pair)



Broda Chair Accessories

- BCM1 seat and back cover
- BCM2L bolster cover left
- BCM2R bolster cover right
- BCM3 neck pillow cover
- BCM4 single foot rest cover
- BCM5 split foot rest cover
- BCM6 armrest covers (pair)



Wheel Chair Pad*

MC1 - Contour Shape 21" by 18"
MC2 - Standard Shape 17" by 15"



Long Split Boot
Adjusts to suit the leg & foot

LSBS (small) LSBL (large)
LSBM (med) LSBXL (X-large)



Short Split Boot
For Hammer Toe & Diabetic ulcers

SSBS (small) SSBL (large)
SSBM (med) SSBXL (X-large)



Medical Sheepskin Heel Dressing

MSHD - 3 per pack
Heel Dressings to be used in conjunction with Heel and Elbow Protectors.



SheepZorb Disposable Dressings (4" x 4")

SZ10 - box of 10
SZ25 - bag of 25



Regular Long Boot
(easy access - 4 sizes)

LBMS (small) LBML (large)
LBMM (med) LBMXL (x-large)



Regular Short Boot
(easy access - 4 sizes)

SBMS (small) SBML (large)
SBMM (med) SBMXL (x-large)



Coccyx Pads (5" by 5")

CP1-Pack of 6
CP2-Pack of 12



Reusable Coccyx Panties and Pads

RCP-2 (small - medium)
RCP-3 (medium -large)



Wrap Around Slipper
with Rubber Sole

EPSMHS (small) EPSMHL (large)
EPSMHM (med) EPSMHXL (X-large)



Pressure Care Slipper
with Rubber Sole

PCBMHS (small) PCBMHL (large)
PCBMHM (med) PCBMHXL (X-large)



Heel Protectors (pair)

UNIHPM



Elbow Protector (1 each)

EPSM



Medical Slippers

MSCR - Medical Sheepskin Cuff with Rubber Sole
MSR - Medical Sheepskin with Rubber Sole



Domestic Slippers

DSCR - Domestic Sheepskin Cuff with Rubber Sole
DSR - Domestic Sheepskin with Rubber Sole



Medical Mules (4 sizes)

SSMS (small) SSML (large)
SSMM (med) SSMXL (X-large)



Palm Protectors (pair)

PPM - Regular
PPMF - with Finger Separators

Pressure relief, moisture absorption and shear reduction you can rely on.

The most naturally wonderful, comfortable and effective product for the prevention and treatment of wounds that health care has ever seen. Once you explore and feel the texture of Australian Medical Sheepskin, you will understand the amazing benefits to both patient and caregiver alike. This product was specifically designed to save health care millions of dollars by not only treating but preventing bed sores, while at the same time giving comfort to long-term and palliative care patients.



The billiard balls show the superior support capacity of Australian Medical Sheepskin as compared to an imitation.

Look for the Leather back...

Only real sheepskin has a leather back. Manufactured fabrics have a knit back that may or may not be coated with a rubber or plasticized material.



Bed Wrap

As a pillow cover, for shoulders, coccyx, or heels and ankles.

BWM1