





the scoop & incremental household measures



HOUSEHOLD MEASURES IN GRAMS

If you do not have the Gelmix scoop, or need to mix smaller increments, you can use a gram scale or one these common household measures.

The Gelmix jar contains a half serving scoop, equivalent to a level 1/2 teaspoon.

GELMIX SCOOP = 1.2 GRAMS

1/2 LEVEL TSP = 1.2 GRAMS



1/4 LEVEL TSP = 0.6 GRAMS

1/8 LEVEL TSP = 0.3 GRAMS



1/16 LEVEL TSP = 0.15 GRAMS

DIGITAL GRAM SCALE

A scale that measures down to 0.01g is the most accurate way to measure small increments of Gelmix. Just "TARE" the scale to bring it back to 0.00 after placing your container on the scale.

