

Performance Data for the Aquasana Drinking Water Filter Models AQ-4000/4600/4601				
Replacement	Operating pressure range	Rated capacity	Operating temp range	Rated flow
AQ-4035	20-80 psi (1.40-5.624 kg/cm²)	450 gallons	40-90° F (4.44-32.2° C)	o.5 gpm
Manufactured by: Aguasana, Inc. 6310 Midway Road · Haltom City, Texas 76117 · 866.662.6885				

Organic chemicals included by surrogate testing Drinking water

regulatory level (MCL/MAC) mg/L

0.002

0.003

0.04

0.005

0.1

0.07

0.6

0.075

0.005

0.007

0.07

0.005

0.007

0.002

0.00005

0.0004

0.05

0.04

0.001

0.004

0.005

0.05

0.07

0.2

0.005

0.005

0.080

10

0.1

0.0002

0.7

0.0002

Fffluent/

Filtered

0.001

0.003

0.001

0.0018

0.001

0.0002

0.0017

0.001

0.001

0.0048

0.0005

0.001

0.001

0.001

0.0002

0.00059 99%

0.00002 >99%

0.0005

0.0006

0.0002

0.0003

0.0001

0.0003

0.0002

0.001

0.000002 >99%

0.00001

0.0001

0.001

0.004

0.0005

0.001

0.001

0.001

0.0016

0.001

0.0005

0.0046 95%

0.0005

0.0010

Effluent/

Filtered

0.015

0.001

0.00001 >99%

0.001

0.00002

Percent

>98%

>97%

>99%

>99%

98%

>99%

99%

98%

>99%

>99%

>98%

95%

>99%

>99%

>99%

>99%

98%

98%

98%

98%

99%

96%

98%

>98%

>99%

>99%

>99%

>97%

>99%

>99%

>99%

>99%

99%

>98%

>99%

>99%

>99%

95%

>99%

Percent Reduction

Reduction

Influent/

Unfiltered

0.050

0.100

0.081

0.190

0.078

0.077

0.015

0.110

0.052

0.080

0.040

0.088

0.083

0.170

0.086

0.080

0.079

0.170

0.053

0.088

0.044

0.022

0.024

0.0006

0.0072

0.0082

0.025

0.044

0.060

0.055

0.050

0.096

0.120

0.150

0.081

0.081

0.078

0.270

0.042

0.160

0.084

0.150

0.180

0.300

0.070

Influent/

Unfiltered

0.015

Testing Performed under NSF/ANSI Standards 42 and 53 and in accordance with the California Department of Health Services Drinking Water Treatment Device Program. This system has been tested according to NSF/ANSI 42, 53 & 401 for reduction of the substances listed below. The concentration of the indicated substances in water entering the system was reduced to a concentration less than or equal to the permissible limit for water leaving the system, as specified in NSF/ANSI 42,53 & 401.

VOCs

chloroform) alachlor

atrazine

benzene

carbofuran

chlorobenzene

chloropicrin

2,4-D

(DBCP) o-dichlorobenzene

dinoseh

endrin

ethylbenzene

carbon tetrachloride

dibromochloropropane

p-dichlorobenzene

1,2-dichloroethane

1.1-dichloroethylene

1,2-dichloropropane

cis-1,2-dichloroethylene

trans-1,2-dichloroethylene

cis-1,3-dichloropropylene

ethylene dibromide (EDB)

haloacetonitriles (HAN) Bromochloroacetontrile

Dibromoacetontrile

Dichloroacetontrile

Trichloroacetontrile

heptachlor epoxide

hexachlorobutadiene

lindane

simazine

styrene

toluene

methoxychlor

pentachlorophenol

1,1,2,2-tetrachloroethane

tetrachloroethylene

tribromoacetic acid

1.2.4-trichlorobenzene

1,1,1-trichloroethane

1.1.2-trichloroethane

Trihalomethanes (THMs)

Bromodichloromethane (THM)

Chlorodibromomethane

trichloroethylene

Bromoform (THM)

Chloroform (THM)

Xylenes (total)

2,4,5-TP (silvex)

1,1,1-trichloro-2-propanone

heptachlor (H-34, Heptox)

hexachlorocyclopentadiene

haloketones (HK) 1,1-dichloro-2-propanone

(by surrogate testing using

NSF/ANSI 42	Influent Challenge Concentration	Reduction requirement	Overall % Reduction	Results
Chlorine Reduction, Free Available	2.0 ± 10% mg/L	50%	>97-4%	Pass
Chloramine Reduction, Free Available	3.0 ± 10% mg/L	.5 mg/L	>97-4%	Pass
Particulate Reduction Class I (0.5 μ to < 1 μ)	at least 10,000 particles/mL	85%	>99.9%	Pass

NSF/ANSI 53	Influent Challenge Concentration	Reduction requirement	Overall % Reduction	Results
Cyst Live Cryptosporidium & Giardia	min 50,000/L	99.95%	>99.99%	Pass
Mercury Reduction pH 8.5	o.oo6 ± 10% mg/L	<2 ug/L	>95.8%	Pass
Mercury Reduction pH 6.5	o.oo6 ± 10% mg/L	<2 ug/L	>96.5%	Pass
Lead Reduction pH 6.5	.15 ± 10% mg/L	<10 ug/L	>99.3%	Pass
Lead Reduction pH 8.5	.15 ± 10% mg/L	<10 ug/L	>99.3%	Pass
MTBE Reduction	.15 ± 20% mg/L	<5 ug/L	86.6%	Pass
Turbidity	11 ± 1 NTU	<0.5 NTU	99.0%	Pass
VOC Surrogate Test	300 ± 30 μg/L	95%	99.4%	Pass
Asbestos Reduction	107 to 108 fibers/L; fibers greater than 10 µm in length	99%	>99%	Pass

NSF/ANSI 401	Maximum Concentration	Minimum Reduction	Overall % Reduction	Results
Phenytoin	30 ng/L	95.50%	95.6%	Pass
Ibuprofen	6o ng/L	95.3%	95.4%	Pass
Naproxen	20 ng/L	96.3%	96.4%	Pass
Estrone	20 ng/L	96.30%	96.5%	Pass
Bisphenol A	300 ng/L	98.80%	98.9%	Pass
Nonyl phenol	200 ng/L	97.50%	97.5%	Pass



System tested and certified by NSF International against NSF/ANSI Standard 42 and 53 for the reduction of the claims specified on the Performance Data Sheet and at www.nsf.org.



System tested and certified by WQA against NSF/ ANSI Standards 42 and 53 for the reduction of the claims specified on the Performance Data Sheet.



Filter is only to be used with cold water.



Filter usage must comply with all state and local laws.



Testing was performed under standard laboratory conditions, actual performancé may vary.

- · All contaminants reduced by this filter are listed.
- · Not all contaminants listed may be present in your water.
- · Filter does not remove all contaminants that may be present in tap water.



Systems certified for cyst reduction may be used on disinfected waters that may contain filterable cysts.



See owner's manual for general installation conditions and needs plus manufacturer's limited warranty.



(THM)

Do not use with water that is microbiologically unsafe or of unknown water quality without adequate disinfection before or after the system.