




Cytek  BIOSCIENCES

Custom Aria Configurations

Custom Aria Configurations

Expand the capabilities of your FACSAria™ by adding lasers and detectors to your trusted Aria cell sorter. Add high power lasers and highly efficient PMTs and update to the latest version of DIVA—while saving money.

Excellent instrument support and performance verification


- Instrument training
- Experienced Service team with guaranteed response times
- CST and QbSure® validation software



Cytek 4-PMT box mounted on Aria door

| Wavelength (nm) | Power (mW) |
|-----------------|------------------|
| 355 | 20,40,60,100 |
| 375 | 16,50,60 |
| 405 | 50,100 |
| 445 | 75 |
| 488 | 20,35,60,100,150 |
| 532 | 50,80,100,150 |
| 552 | 50,80,100,150 |
| 561 | 50,80,100,150 |
| 594 | 20,60,100,200 |
| 640 | 25,40,100 |

Available lasers



Lasers on front optical plate

About Cytek

We are dedicated to providing the flow cytometry tools scientists need at prices they can afford. Cytek has training and technical support available from many of the leading engineers in Flow Cytometry analyzer and sorter instruments.