



Cytek  BIOSCIENCES

Custom Canto Configurations

Custom Canto Configurations


Enhance your FACSCanto™ for less by adding lasers and detectors, Add customized lasers and highly efficient PMTs for a flexible configuration, and get the latest version of DIVA.

Wavelength (nm)	Power (mW)
405	50,100
445	45
488	20,40,60,80,120
532	20,40,60,80,100,120
552	40,60,80,100,120
561	40,60,80,100,120
594	40
640	17,75


Available lasers

Excellent instrument support and performance verification

- Instrument training
- Experienced Service team with guaranteed response times
- CST and QbSure® validation software



Bead data from FACSCanto 2 Violet laser upgrade



Lasers on Canto front optical plate

About Cytek

We are dedicated to providing the flow cytometry tools scientists need at prices they can afford. Cytek has training and technical support available from many of the leading engineers in Flow Cytometry analyzer and sorter instruments.